
**Prevena**<sup>TM</sup>  
Incision Therapy

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
Publikation	Wund-/ Operationstyp	Klinischer Evidenzgrad*
Newman JM, Siqueira MBP, Klika AK, Molloy RM, Barsoum WK, Higuera CA. Use of Closed Incisional Negative Pressure Wound Therapy After Revision Total Hip and Knee Arthroplasty in Patients at High Risk for Infection: A Prospective, Randomized Clinical Trial. <i>Journal of Arthroplasty</i> . 2018 Nov 17. [Epub Ahead of Print]	Hüft- und Knie-Totalendoprothese	1b
Crist BD, Oladeji LO, Khazzam M, Della Rocca GJ, Murtha YM, Stannard JP. Role of acute negative pressure wound therapy over primarily closed surgical incisions in acetabular fracture ORIF: A prospective randomized trial. <i>Injury</i> . 2017 Apr 27;pii: S0020-1383(17)30283-8.	Acetabulum-frakturen	1b
Pauser J, Nordmeyer M, Biber R, Jantsch J, Kopschina C, Bail HJ, Brem MH. Incisional negative pressure wound therapy after hemiarthroplasty for femoral neck fractures - reduction of wound complications. <i>International Wound Journal</i> . 2014;13(5):663-667.	Hemiprothese bei Oberschenkelhalsfrakturen	1b
Manoharan V, Grant A, Harris A, Hazratwala K, Wilkinson M, McEwen P. Closed Incision Negative Pressure Wound Therapy vs Conventional Dry Dressings After Primary Knee Arthroplasty: A Randomized Controlled Study. <i>J Arthroplasty</i> . 2016 Apr 28. pii: S0883-5403(16)30083-3.	Knie-Endoprothese	1b
Howell RD, Hadley S, Strauss E, Pelham FR. Blister formation with negative pressure dressings after total knee replacement. <i>Current Orthopaedic Practice</i> . 2011 Mar;22(2):176-179.	Knie-Endoprothese	1b
Stannard JP, Robinson JT, Anderson ER, McGwin G Jr, Volgas DA, Alonso JE. Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma. <i>Journal of Trauma</i> . 2006 Jun;60(6):1301-6.	Frakturen der unteren Extremitäten	1b
Stannard JP, Volgas DA, McGwin G, Stewart RL, Obremsky W, Moore T, Anglen JO. Incisional negative pressure wound therapy after high-risk lower extremity fractures. <i>Journal of Orthopedic Trauma</i> . 2012 Jan;26(1):37-42.	Frakturen der unteren Extremitäten	1b
Stannard JP, Volgas DA, Stewart R, McGwin G Jr, Alonso JE. Negative pressure wound therapy after severe open fractures: a prospective randomized study. <i>Journal of Orthopedic Trauma</i> . 2009 Sep;23(8):552-7.	Frakturen der unteren Extremitäten	1b
Pachowsky M, Gusinde J, Klein A, Lehl S, Schulz-Drost S, Schlechtweg P, Pauser J, Gelse K, Brem MH. Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty. <i>International Orthopaedics</i> . 2012 Apr;36(4):719-22.	Hüft-Totalendoprothese	1b
Redfern RE, Cameron-Ruetz C, O'Drobinak S, Chen J, Beer KJ. Closed incision negative pressure therapy effectson postoperative infection and surgical site complication after total hip and knee arthroplasty. <i>J Arthroplasty</i> 2017 Nov;32(11):3333-3339.	Endoprothese Knie und Hüfte	2
Reddix RN Jr, Leng XI, Woodall J, Jackson B, Dedmond B, Webb LX. The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery. <i>Journal of Surgical Orthopaedic Advances</i> . 2010 Jun;19(2):91-7.	Hüft-Endoprothese	3
Cooper HJ, Roc GC, Bas MA, Berliner ZP, Hepinstall MS, Rodriguez JA, Weiner LS. Closed incision negative pressure therapy decreases complications after periprosthetic fracture surgery around the hip and knee. <i>Injury</i> . 2018 Feb;49(2):386-391. doi: 10.1016/j.injury.2017.11.010. Epub 2017 Nov 14.	Periprosthetische Fraktur-Chirurgie	3
Cooper HJ, Bas MA. Closed-Incision Negative-Pressure Therapy Versus Antimicrobial Dressings After Revision Hip and Knee Surgery: A Comparative Study. <i>J Arthroplasty</i> . 2016 May;31(5):1047-52.	Revisionseingriff Knie und Hüfte	3
Anatone AJ, Shah RP, Jennings EL, Geller JA, Cooper J. A risk-stratification algorithm to reduce superficial surgical site complications in primary hip and knee arthroplasty. <i>Arthroplasty Today</i> . 2018;4(4):493-498. doi:10.1016/j.artd.2018.09.004.	Endoprothese Knie und Hüfte	3
Curley AJ, Terhune EB, Velott AT, Argintar EH. Outcomes of Prophylactic Negative Pressure Wound Therapy in Knee Arthroplasty. <i>Orthopaedics</i> . 2018;41(6):e837-e840. doi:10.3928/01477447-20181010-02.	Knie-Endoprothese	3

Publikation	Wund-/ Operationstyp	Klinischer Evidenzgrad*
Reddix RN, Tyler HK, Kulp B, Webb LX. Incisional vacuum-assisted wound closure in morbidly obese patients undergoing acetabular fracture surgery. <i>The American Journal of Orthopedics</i> . 2009 Sep;38(9):32-5.	Acetabulum-frakturen	4
Hansen E, Durinka JB, Costanzo JA, Austin MS, Deirmengian GK. Negative pressure wound therapy is associated with resolution of incisional drainage in most wounds after hip arthroplasty. <i>Clinical Orthopaedics and Related Research</i> . 2013 Oct;471(10):3230-6.	Hüft-Endoprothese	4
Stannard JP, Atkins BZ, O-Malley D, Singh H, Bernstein B, Fahey M, Masden D, Attinger CE. Use of negative pressure therapy on closed surgical incisions: A case series. <i>Ostomy Wound Management</i> . 2009 Aug;55(8):58-66.	Frakturen der unteren Extremitäten	4
Gomoll AH, Lin A, Harris MB. Incisional vacuum-assisted closure therapy. <i>Journal of Orthopaedic Trauma</i> . 2006 Nov-Dec;20(10):705-9.	Orthopädisches Trauma	4
Stannard JP, Gabriel A, Lehner B. Use of Negative Pressure Wound Therapy Over Clean, Closed Surgical Incisions. <i>International Wound Journal</i> . 2012;9:32-39.	Orthopädisches Trauma	4
Berkowitz MJ. Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures. <i>Techniques in Foot and Ankle Surgery</i> . 2013 Dec;12(4):172-174.	Kalkaneusfrakturen	5
Brem MH, Bail HJ, Biber R. Value of Incisional Negative Pressure Wound Therapy in Orthopedic Surgery. <i>International Wound Journal</i> . 2014 Jun;11(Suppl 1):3-5.	Gemischt	5
Suleiman LI, Mesko DR, Nam D. Intraoperative Considerations for Treatment/Prevention of Prosthetic Joint Infection. <i>Current Reviews in Musculoskeletal Medicine</i> . 2018:1-8.	Endoprothese Knie und Hüfte	5
Chotanaphuti T, Courtney PM, Fram B, Kleef N.J., Kim TK, Kuo FC, Lustig S, Moojen DJ, Nijhof M, Oliashirazi A, Poolman R, Purtill JJ, Rapisarda A, Rivero-Boschert S, Veltman ES. Hip and Knee Section, Treatment, Algorithm: Proceedings of International Consensus on Orthopedic Infections. <i>The Journal of Arthroplasty</i> . 34(2S):S393-S397. doi: 10.1016/j.arth.2018.09.024.	Endoprothese Knie und Hüfte	5
DeCarbo WT, Hyer CF. Negative-Pressure Wound Therapy Applied to High-Risk Surgical Incisions. <i>Journal of Foot and Ankle Surgery</i> . 2010 Mai;49(3):299-300.	Orthopädisches Trauma	5
Nam D, Sershon RA, Levine BR, Della Valle CJ. The Use of Closed Incision Negative-Pressure Wound Therapy in Orthopaedic Surgery. <i>J Am Acad Orthop Surg</i> . 2018:1-8. doi: 10.5435/JAAOS-D-17-00054.	Orthopädische Chirurgie	5
Al-Hourai RK, Aalirezaie A, Adib F, Anoushiravani A, Bhashyam A, Binlaksar R, Blevins K, Bonanzinga T, Chih-Kuo F, Cordova M, Deirmengian GK, Fillingham Y, Frenkel T, Gomez J, Gundtoft P, Harris MA, Harris M, Heller S, Jennings JA, Jimenez-Garrido C, Karam JA, Khlopas A, Klement MR, Komnos G, Krebs V, Lachiewicz P, Miller AO, Mont MA, Montanez E, Romero CA, Schwarzkopf R, Shaffer A, Sharkey PF, Smith BM, Sodhi N, Thienpont E, Villanueva AO, Yazdi H. General Assembly, Prevention, Wound Management: Proceedings of International Consensus on Orthopedic Infections. <i>The Journal of Arthroplasty</i> . 2019;34(2):S157-S168. doi:10.1016/j.arth.2018.09.066.	Orthopädische Infektionen	5
Agarwal A. Management of Closed Incisions Using Negative-Pressure Wound Therapy in Orthopedic Surgery. <i>Plastic and reconstructive surgery</i> . 2019;143(1S Management of Surgical Incisions Utilizing Closed Incision Negative Pressure Therapy):21S-26S.	Orthopädische Unfallchirurgie	5

## Referenzen

- Public Health England. Surveillance of surgical site infections in NHS hospitals in England. April 2018 bis März 2019. Veröffentlicht im Dezember 2019.
- Jenks, P.J. Clinical and economic burden of surgical site infection (SSI) and predicted financial consequences of elimination of SSI from an English hospital, Volume 86 (2014), Issue 1, pg 24–33.
- Merkow R. Underlying reasons associate with hospital readmission following surgery in the US. 2015;313(5):483–95.
- M. Nobile, P. Navone, A. Orzella, *et al.* Developing a model for analysis the extra costs associated with surgical site infections (SSIs): an orthopaedic and traumatological study run by the Gaetano Pini Orthopaedic Institute, 4 (2015), p. P68.
- Wilkes RP, Kilpadi DV, Zhao Y. Closed Incision Management With Negative Pressure Wound Therapy (CIM): Biomechanics. 2012;19(1):67–75.
- Kilpadi DV, Cunningham MR. Evaluation of Closed Incision Management with Negative Pressure Wound Therapy (CIM): Hematoma/Seroma and Involvement of the Lymphatic System. 2011;19:588–596.
- Glaser DA, Farnsworth CL, Varley ES. Negative pressure therapy for closed spine incisions: A pilot study. 2012;24(11):308–316.
- Federal Drug Administration. De Novo Classification Request for PREVENA 125 and PREVENA PLUS 125 Therapy Units. De Novo Summary (DEN180013), 2019. [https://www.accessdata.fda.gov/cdrh\\_docs/reviews/DEN180013.pdf](https://www.accessdata.fda.gov/cdrh_docs/reviews/DEN180013.pdf)
- Newman JM, Siqueira MBP, Klika AK, Molloy RM, Barsoum WK, Higuera CA. Use of Closed Incisional Negative Pressure Wound Therapy After Revision Total Hip and Knee Arthroplasty in Patients at High Risk for Infection: A Prospective, Randomized Clinical Trial. *J Arthroplasty*. 2019;34(3):554–559.
- Cooper HJ, Bas MA. Closed-Incision Negative-Pressure Therapy Versus Antimicrobial Dressings After Revision Hip and Knee Surgery: A Comparative Study. *J Arthroplasty*. 2016;31(5):1047–1052.
- Anatone AJ, Shah RP, Jennings EL, Geller JA, Cooper J. A risk-stratification algorithm to reduce superficial surgical site complications in primary hip and knee arthroplasty. *Arthroplasty Today*. 2018;4(4):493-498. doi:10.1016/j.artd.2018.09.004.

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© 3M 2022. 3M, ActiV.A.C., Arthro-Form, Prevena, Prevena Restor und V.A.C. sind eingetragene Marken oder Warenzeichen der 3M Company. Alle Rechte vorbehalten.
**3M Deutschland GmbH**
**3M Österreich GmbH**
**3M (Schweiz) GmbH**