Proaktives Risikomanagement mit der Unterdrucktherapie für geschlossene Inzisionen (Closed Incision Negative Pressure Therapy, ciNPT)

Ein klinisches ciNPT-Kompendium von 3M



Brief des Herausgebers

Liebe Kolleginnen und Kollegen,

Wir freuen uns sehr, Ihnen dieses Dokument vorstellen zu können. Es enthält evidenzbasierte Leitlinien zum Einsatz der 3M™ Prevena™ Therapie für das proaktive Risikomanagement bei einer Vielzahl spezifischer klinischer Situationen. Die Unterdruckwundtherapie zählt zu den wichtigsten Innovationen bei der Wundversorgung in den letzten 30 Jahren. In jüngster Zeit hat sich der Einsatz der Unterdrucktherapie für geschlossene chirurgische Inzisionen als verlässliche Methode etabliert, um das Risiko postoperativer Komplikationen bei Patienten, bei denen ein solches Risiko besteht, zu reduzieren.

Bei der Prevena Therapie wird eine Unterdrucktherapie auf geschlossene Inzisionen und das umgebende Weichgewebe angewandt, um das Behandlungsergebnis zu optimieren und Komplikationen zu reduzieren. Eine der häufigsten Anfragen, die ich von Chirurgen erhalte, sind spezifische Empfehlungen, wann genau sie die Prevena Therapie bei ihren Patienten anwenden sollten. Mit mehr als 200 von Experten geprüften Veröffentlichungen zur Prevena Therapie liegen nun Erkenntnisse in ausreichendem Umfang vor, um evidenzbasierte Leitlinien zu erstellen, die Chirurgen bei ihrer Entscheidungsfindung helfen können. Natürlich sind diese Leitliniendokumente nicht als Ersatz für Entscheidungen auf der Grundlage von klinischem Fachwissen gedacht. Sie sollen den Chirurgen nur als zusätzliche Gesichtspunkte dienen, die auf der aktuellen veröffentlichten Literatur basieren.

Wir sind der Auffassung, dass die konsequente Anwendung der Prevena Therapie bei geeigneten Patienten für proaktives Risikomanagement zu besseren Behandlungsergebnissen beitragen sowie das Risiko von Komplikationen und die Gesamtbehandlungskosten mindern kann. Wir hoffen, dass dieses Dokument für Sie von Nutzen ist.

Mit freundlichen Grüßen,

Ron Silverman, MD

Senior Vice President of Clinical Affairs und Chief Medical Officer, 3M Health Care Business Group



Ronald P. Silverman, MD, FACS

Senior Vice President, Global Medical and Clinical Affairs, Chief Medical Officer, 3M Health Care Business Group Clinical Associate Professor of Surgery, University of Maryland Adjunct Associate Professor of Plastic Surgery, Johns Hopkins

Dr. Ron Silverman ist Senior
Vice President, Global Medical and
Clinical Affairs und Chief Medical
Officer für die 3M Health Care
Business Group. Er ist außerdem
praktizierender plastischer und
rekonstruktiver Chirurg am University
of Maryland Medical Center.
Vor seiner Rolle als CMO war
Dr. Silverman Chief of Plastic Surgery
an der University of Maryland und
folgte akademischen Rufen an die
University of Maryland und die Johns
Hopkins Schools of Medicine.

Ron schloss seine Facharztausbildung für plastische Chirurgie am Harvard/ Massachusetts General Hospital ab und machte seinen Abschluss mit Cum Laude an der University of Maryland School of Medicine. Er wurde vom American Board of Plastic Surgery zertifiziert und ist Mitglied verschiedener Berufsverbände, einschließlich des American College of Surgeons, der American Society of Plastic Surgeons, des Plastic Surgery Research Council und der American Association of Plastic Surgeons. Er erhielt den Wayne W. Babcock Award für herausragende Leistungen in der Chirurgie und ist Mitglied der Alpha Omega Alpha Honor Society.

PRM Proaktives Risikomanagement (PRM) mit der 3M™ Prevena™ Therapie

Inhaltsverzeichnis

04	Übersicht über das proaktive Risikomanagement (PRM)	
05	Wirkmechanismus der Prevena Therapie	
06	ciNPT-Metaanalyse 2023: verschiedene Fachgebiete	
80	Identifizieren von Patienten und Verfahren, die von der Prevena Therapie profitieren können	
	beachten Sie, dass die folgenden Kapitel in der Originalsprache Englisch sst sind.	
09	PRM in der Orthopädie	
27	PRM in der plastischen Chirurgie	
44	PRM in der Gefäßchirurgie	
57	PRM in der Herz-Thorax-Chirurgie	
67	PRM in der Allgemeinchirurgie	
73	Produktübersicht	
75	Kontaktinformationen	

Schützen Sie Ihre Patienten mit der 3M™ Prevena™ Therapie.

Einführung des proaktiven Risikomanagements (PRM)

Die Prevena Therapie kann für alle Patienten von Nutzen sein – insbesondere bei Hochrisikopatienten kann sie zu einer Risikominderung bezüglich postoperativer Wundinfektionen* beitragen und die Realisierung von Kosteneinsparungen ermöglichen. Durch die Einführung von PRM können Sie von einer Stratifizierung Ihrer Verfahren und des Patientenrisikos profitieren und so Ihre Hochrisikopatienten besser schützen.

Komplikationen an der Eingriffsstelle (Surgical Site Complications, SSC) sind nicht nur kostspielig, sie können auch die Genesung der Patienten beeinträchtigen.



Postoperative Wundinfektionen (Surgical Site Infections, SSIs) treten bei 2 - 5 % aller Patienten auf.1



Bei Patienten, bei denen eine SSI auftritt, ist die Wahrscheinlichkeit einer erneuten Aufnahme etwa fünfmal höher.2



Durchschnittlich über 5.000 € zusätzliche Kosten pro SSI.3

- * Die Wirksamkeit der Prevena Therapie bei der Reduzierung der Inzidenz von postoperativen Wundinfektionen und Seromen wurde nicht bei allen chirurgischen Eingriffen und Patientengruppen nachgewiesen. Die vollständigen Indikationen zur Anwendung und Einschränkungen finden Sie auf hcbgregulatory.3M.com
- † Die Berechnungen wurden basierend auf der relativen Patientengruppen-Inzidenz, die in dieser Studie berichtet wurde, abgeleitet.
- ‡ Statistisch signifikant (p≤0,05).

Im Rahmen der PROMISES-Studie (Post-market, Randomized, Open-Label, Multicenter Study to evaluate Effectiveness, zu deutsch: nach der Markteinführung durchgeführte offene, multizentrische, randomisierte Studie zur Bewertung der Wirksamkeit) wurde die Wirksamkeit der Unterdrucktherapie bei geschlossenen Inzisionen im Vergleich zu silberhaltigen Wundauflagen in Bezug auf die Reduzierung postoperativer Komplikationen bei Hochrisikopatienten gemessen, die sich der Revision einer Knieendoprothese unterzogen hatten.⁴

Die 3M™ Prevena™ Therapie reduziert nachweislich das Risiko von SSCs sowie die Gesamtbehandlungskosten. 4,5

Die Prevena Therapie zeigte positive Behandlungsergebnisse über mehrere Fachgebiete hinweg, einschließlich der plastischen Chirurgie, der Gefäß-, Herz-Thorax-, oder Wirbelsäulenchirurgie sowie der Orthopädie und Allgemeinchirurgie.⁶ Daten aus einer multizentrischen randomisiert-kontrollierten Studie und gesundheitsökonomischen Analyse zeigten, dass die Prevena Therapie das Risiko von postoperativen Wundinfektionen über 90 Tage⁴ sowie Wiederaufnahmen⁴ und die Kosten für die Versorgung der Eingriffsstelle⁵ im Vergleich zu silberhaltigen Wundauflagen signifikant reduziert.

Reduzierung von SSCs^{†4}

3.4 % (5/147) Prevena Therapie vs. 14.3 % (21/47) derzeitiger Behandlungsstandard $(p=0,0013)^{t}$

Reduzierung der Wiederaufnahmeraten^{†4}

3,4 % (5/147) Prevena Therapie vs. 10,2 % (15/47) derzeitiger Behandlungsstandard (p=0.0208)


Durschnittlich weniger Verbandswechsel^{†4}

1,1 ± 0,29 Prevena Therapie vs. 1,3 ± 0,96 derzeitiger Behandlungsstandard $(p=0,0003)^{1}$

Literatur

1. Anderson, DJ, et al. Strategies to Prevent Surgical Site Infections in Acute Care Hospitals: 2014 Update. Infect Contol Hosp Epidemiol. 2014;35(6):605-627. doi: 10.1086/676022. 2. Canadian Surgical Site Infection Prevention Audit Month Report. Heruntergeladen von http://www.patientsafetyinstitute.ca/en/ toolsResources/ Pages/SSI-Audit-Recap-Report-2016-12.aspx 3. Jenks PJ, Laurent M, McQuarry S, Watkins R. Clinical and economic burden of surgical site infection (SSI) and predicted financial consequences of elimination of SSI from an English hospital. J Hosp Infect. 2014 Jan;86(1):24-33. 4. Higuera-Rueda CA, Emara AK, Nieves-Malloure Y, Klika AK, Cooper HJ, Cross MB, Guild GN, Nam D, Nett MP, Scuderi GR, Cushner FD, Piuzzi NS, Silverman RP. The Effectiveness of Closed-Incision Negative-Pressure Therapy Versus Silver-Impregnated Dressings in Mitigating Surgical Site Complications in High-Risk Patients After Revision Knee Arthroplasty: The PROMISES Randomised Controlled Trial. J Arthroplasty. Jul. 2021;36(7S):S295-S302.e14. FREI ZUGÄNGLICH Beachten Sie, dass die Dauer der Therapie außerhalb des in der Gebrauchsanleitung empfohlenen Zeitraums liegen kann. 5. Cooper HJ, Bongards C, Silverman RP. Cost-Effectiveness of Closed Incision Negative Pressure Therapy [PREVENA] for Surgical Site Management After Revision Total Knee Arthroplasty: Secondary Analysis of a Randomized Clinical Trial. Journal of Arthroplasty. Aug. 2022;37(8S):S790-S795. FREI ZUGÄNGLICH 6. Cooper HJ, Roc GC, Bas MA, et al. Injury. 218;49(2):386-391.

Wirkmechanismus der 3M™ Prevena™ Therapie.



Weitere, einzigartige Vorteile der 3M™ Prevena™-Wundauflagen:



Unter Unterdruck kollabiert die Wundauflage aus retikuliertem, offenporigem Schaumstoff in ihren geometrischen Mittelpunkt



Gleichmäßige Verteilung des Unterdrucks



Die Kontaktschicht enthält 0,019 % ionisches Silber



Erhältlich in verschiedenen Größen und Konfigurationen, die sich für verschiedene Anwendungen und Patienten eignen

Die 3M™ Prevena™ Therapie kann dazu beitragen, dass Ärzte Patienten früher in eine Behandlung in häuslicher Umgebung entlassen können:

- ► Tragbare Therapie zur Einmalverwendung für 7 oder 14 Tage[§]
- ► Zum Duschen geeignet[‡]
- ► Akustische und visuelle Alarme
- ► Dedizierte klinische Unterstützung

HINWEIS: Für diese Produkte und Therapien gibt es spezifische Indikationen, Vorsichtsmaßnahmen, Kontraindikationen, Warnhinweise und Sicherheitsinformationen. Bitte kontaktieren Sie vor der Anwendung medizinisches Fachpersonal und lesen Sie die Bedienungsanleitung.

Literatur

Grauhan O, Navasardyan A, Hofmann M, Muller P, Stein J, Hetzer R. Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy. J Thorac Cardiovasc Surg. 2013;145:1387-1392.
 Kilpadi DV, Cunningham MR. Evaluation of closed incision management with negative pressure wound therapy (CIM): Haematoma/seroma and involvement of the lymphatic system. Wound Repair and Regeneration. 2011;19(5):588-96.
 Glasser, et al. Negative pressure therapy for closed spine incisions: a pilot study. Wounds. 2012;24(11):308-16.
 Wilkes RP, Kilpadi DV, Zhao Y, Kazala R, McNulty A. Closed incision management with negative pressure wound therapy (CIM): biomechanics. Surg Innov. 1. März 2012;19(1):67-75.
 Colli A. First experience with a new negative pressure incision management system on surgical incisions after cardiac surgery in high risk patients. Journal of Cardiothoracic Surgery.
 Dezember 2011;6(1):160.

^{*} In einem Kanister.

[†] In Computer-Referenzmodellen.

[‡] Weitere Informationen finden Sie in der Bedienungsanleitung für Prevena Therapie-Patienten und -Anwender

[§] Die maximale Therapiedauer mit der 3M™ Prevena™ Therapieplattform beträgt 7 Tage. Die maximale Therapiedauer mit der 3M™ Prevena Restor™ Therapieplattform beträgt 14 Tage.

Unterdrucktherapie bei geschlossenen Inzisionen im Vergleich zum Behandlungsstandard zur Verringerung von postoperativen Komplikationen: ein systematischer Review und Meta-Analyse.

Cooper, H. John MD, Singh, Devinder P. MD, Gabriel, Allen MD, FACS, Mantyh, Christopher MD, Silverman, Ronald MD, Griffin, Leah MS. Closed Incision Negative Pressure Therapy versus Standard of Care in Reduction of Surgical Site Complications: A Systematic Review and Meta-analysis. *Plastic & Reconstructive Surgery-Global Open* 11(3):p e4722, March 2023.

Hintergrund

- Postoperative Komplikationen (Surgical Site Complications, SSC) wie postoperative Wundinfektionen (Surgical Site Infections, SSI), Dehiszenz, Serome, Hämatome und Hautnekrosen können die Behandlungsergebnisse beeinträchtigen und sich negativ auf die Behandlungskosten auswirken.
- ► Zur Reduzierung des Risikos postoperativer Komplikationen wurden Optionen zur Behandlung des Operationssitus entwickelt, einschließlich der Unterdrucktherapie für geschlossene Inzisionen (Closed Incision Negative Pressure Therapy, ciNPT*).
- ► Die Verwendung von ciNPT wurde mit positiven Behandlungsergebnissen in vielen chirurgischen Fachgebieten assoziiert.¹-6

Studienzweck

Diese systematische Übersicht und Meta-Analyse bewertet den Einfluss der ciNPT auf postoperative und gesundheitsökonomische Ergebnisse anhand der veröffentlichten Studien.


Methoden


- ► Es wurde eine systematische Literaturrecherche über PubMed, EMBASE und QUOSA durchgeführt.
- Bewertet wurden englischsprachige Publikationen, welche die ciNPT zwischen Januar 2005 und August 2021 mit dem Behandlungsstandard bezüglich Wundauflagen vergleichen.
- Extrahiert wurden die Charakteristika der Studienteilnehmer, die chirurgischen Verfahren, die verwendeten Verbände, die Behandlungsdauer, die postoperativen Behandlungsergebnisse und die Folgeergebnisse.

Ergebnisse

- ► Bei der Literaturrecherche wurden 84 Studien für die Analyse identifiziert.
- ► Es wurden bezüglich der ciNPT signifikante Reduktionsraten in Bezug auf SSC festgestellt (p<0,001).
- ► Signifikante Verringerung von SSI (p<0,001), oberflächliche SSI (p<0,001), tiefe SSI (p=0,002), Serome (p=0,002), Dehiszenz (p=0,022) und Hautnekrosen (p=0,001) wurden mit der Verwendung von ciNPT in Verbindung gebracht (p<0,05).
- ► Bei der ciNPT zeigten sich signifikante Reduktionen bezüglich der Wiederaufnahmen und Revisionseingriffen (p<0,05).
- Bei ciNPT-Patienten war der Krankenhausaufenthalt um 0,9 Tage kürzer als bei Patienten, die mit dem Behandlungsstandard behandelt wurden (p<0,001).
- ► Die Unterschiede bei den postoperativen Schmerzwerten und berichteten Mengen an Opioideinsatz waren signifikant zugunsten der Verwendung von ciNPT (p<0,05).
- ► Auch die Anwendung von postoperativen Drainagen und Antibiotika war bei ciNPT-Patienten verringert, jedoch nicht in signifikantem Ausmaß.

^{*3}M™ Prevena™ Incision Management System





Fazit

- ► Bei diesen Metaanalysen war die Anwendung der ciNPT mit einer statistisch signifikanten Reduktion bei der Inzidenz von postoperativen Komplikationen, postoperativen Wundinfektionen, Seromen, Dehiszenzen und Hautnekrosen assoziiert.
- ► Bei ciNPT-Patienten wurde auch eine Reduktion der Wiederaufnahmen, Revisionseingriffen und der Dauer des Krankenhausaufenthalts beobachtet sowie auch eine Verringerung der Schmerzen und des Einsatzes von Opiaten.
- ► Die Studie wird durch die Kombination aus Beobachtungstudien und randomisierten, kontrollierten Studien sowie durch die Vermischung verschiedener chirurgischer Fachgebiete und die Unterschiede bei den Datenangaben bei den einbezogenen Artikeln eingeschränkt.
- ► Es ist zu beachten, dass sich die Daten auf ein bestimmtes im Handel erhältliches ciNPT-System beziehen und aufgrund der Unterschiede zwischen den Medizinprodukten nicht auf andere verfügbare Systeme anwendbar sind.
- ► Chirurgen sollten sämtliche verfügbaren Daten in Betracht ziehen, bevor sie sich entscheiden, ob sie ein bestimmtes ciNPT-Medizinprodukt verwenden wollen.

*HINWEIS: Hämatome erreichten keine Signifikanz, zeigten aber eine Tendenz, die für die Behandlung spricht

Die Wirksamkeit der Prevena Therapie bei der Reduzierung der Inzidenz von postoperativen Wundinfektionen und Seromen wurde nicht bei allen chirurgischen Eingriffen und Patientengruppen nachgewiesen. Die vollständigen Indikationen zur Anwendung und Einschränkungen finden Sie auf hcbgregulatory.3M.com.

Literatur

Cooper HJ, Roc GC, Bas MA, et al. *Injury*. 218;49(2):386-391.
 Ruggieri VG, Olivier ME, Aludaat C, et al. *Heart Surg Forum*. 2019;22(2):E092-E096.
 Licari L, Campanella S, Carolla C, et al. *Cureus*. 2020;12(5):e8283.
 Swift SH, Zimmerman MB, Hardy-Fairbanks AJ. *J Reprod Med*. 2015;60(5-6):211-218.
 Ferrando PM, Ala A, Bussone R, et al. *Plast Reconstr Surg Glob Open*. 2018;6(6):e1732.
 Pleger SP, Nink N, Elzien M, et al. *Int Wound J*. 2018;15(1):75-83.

¹Columbia University Irving Medical Centre, New York, New York; ²University of Miami Health System and Miller School of Medicine, Miami, FL; ³Plastic Surgery, Vancouver, WA; ⁴Duke University School of Medicine, Durham, NC; ⁵University of Maryland School of Medicine, Baltimore, MD; ⁶3M, St Paul, MN

<u>Hier die vollständige Studie lesen</u>



Herausgeber: Plastic and Reconstructive Surgery-Global Open

Titel: Closed Incision Negative Pressure Therapy versus Standard of Care in Reduction of Surgical Site Complications: A Systematic Review and Meta-analysis.

Veröffentlicht: März 2023

Cooper, H. John MD, Singh, Devinder P. MD, Gabriel, Allen MD, FACS, Mantyh, Christopher MD, Silverman, Ronald MD, Griffin, Leah MS. Closed Incision Negative Pressure Therapy versus Standard of Care in Reduction of Surgical Site Complications: A Systematic Review and Meta-analysis.

Plastic & Reconstructive Surgery-Global Open 11(3):p e4722, March 2023. | DOI: 10.1097/GOX.000000000000004722 FREI ZUGÄNGLICH

Patienten und Verfahren, die von der 3M™ Prevena™ Therapie profitieren können.1

Eine interdisziplinäre Gruppe von Chirurgieexperten und Fachleuten für Infektionskrankheiten hat einen Algorithmus entwickelt, um zu entscheiden, wann eine Unterdrucktherapie bei geschlossenen Inzisionen (3M™ Prevena™ Therapie) angewandt werden sollte.

Die Konsensempfehlungen basieren auf:


- ► Einer Literaturprüfung
- Erfahrungen mit der ciNPT
- Bekannte Risikofaktoren bezüglich postoperativer Vorfälle (Surgical Site Occurrences, SSOs)

Erkenntnisse:

- ► Zahlreiche Veröffentlichungen berichten von Risikofaktoren in Bezug auf postoperative Wundinfektionen, wobei zu den häufigsten Adipositas (Body-Mass-Index ≥30 kg/m²), Diabetes mellitus, Nikotinabusus oder eine verlängerte Operationszeit zählen.
- ► Chirurgen wird empfohlen, die individuellen Risikofaktoren und chirurgischen Risiken eines Patienten zu bewerten.

Chirurgen können die Anwendung der ciNPT bei Patienten in Betracht ziehen, bei denen ein hohes Risiko für das Auftreten von SSO besteht oder die sich einer Prozedur unterziehen, die stark risikobehaftet ist oder die eine hohe Morbidität nach sich ziehen könnte, falls es zu einer postoperativen Wundinfektion kommen sollte.

Bewertung der Risikofaktoren bei einer Unterdrucktherapie für geschlossene Inzisionen (ciNPT):



Weitere zu berücksichtigende Faktoren:

Patientenbezogene Risikofaktoren

- ► Diabetes mellitus
- ► ASA-Wert ≥3
- ► Fortgeschrittenes Alter
- Adipositas
- ► Aktiver Tabakkonsum
- ► Hypoalbuminämie
- Anwendung von Kortikosteroiden
- Aktiver Alkoholkonsum
- Geschlecht
- Chronische Niereninsuffizienz
- Chronisch obstruktive Lungenerkrankung

Allgemeine inzisionsbezogene Risikofaktoren

- ► Inzision unter hoher Spannung
- ► Wiederholte Inzisionen
- ► Umfassende Unterminierung
- ► Traumatisiertes Weichteilgewebe ► Notoperation
- ▶ Ödem
- ► Kontamination
- Verlängerte Operationszeit
- ► Postoperative Bestrahlung
- ► Mechanisch unvorteilhafte Eingriffsstelle

Verfahrens-/operationsbezogene Risikofaktoren

Allgemeinchirurgie

- ► Offene Operation
- ► Kolorektale Operation
- ► Urologische Operation
- ► Behandlung eines Narbenbruchs

Plastische Chirurgie

- ► Bauchdeckenrekonstruktion
- ► Brustrekonstruktion
- ► Große Weichteildefekte
- ► Gynäkologische Operation ► Kontaminierungsrisiko

Orthopädie

- ► Offene Reposition und interne Fixierung von Frakturen
- ▶ Fasziotomie
- Amoutation oberhalb/ unterhalb des Knies

Gefäßchirurgie

- ► Amputation oberhalb/ unterhalb des Knies
- ► Implantierung synthetischer Transplantate

Herzchirurgie

► Sternotomie

Literatur

1. Willy C, Agarwal A, Andersen CA, et al. Closed incision negative pressure therapy: international multidisciplinary consensus recommendations. Int Wound J. 2017;14(2):385-398. doi:10.1111/iwj.12612.



Hier die vollständige Studie lesen





PROMISES study Data PRM in orthopaedic surgery

PROMISES study data suggests 3M™ Prevena™ Therapy can help advance the standard of care.

Higuera-Rueda CA, Emara AK, Nieves-Malloure Y, Klika AK, Cooper HJ, Cross MB, Guild GN, Nam D, Nett MP, Scuderi GR, Cushner FD, Piuzzi NS, Silverman RP. The Effectiveness of Closed-Incision Negative-Pressure Therapy Versus Silver-Impregnated Dressings in Mitigating Surgical Site Complications in High-Risk Patients After Revision Knee Arthroplasty: The PROMISES Randomised Controlled Trial. J Arthroplasty. 2021 Jul;36(7S):S295-S302.e14.

The PROMISES (Post-market, Randomised, Open-Label, Multicentre study to evaluate Effectiveness) Trial.

Data from a multicentre randomised controlled trial showed that Prevena Therapy significantly reduced the risk of 90-day surgical site complications (SSCs) and post-op readmissions vs. silver-impregnated dressings.

Study design

Post-market, randomised, open-label, multicentre study (United States).

Study purpose

Evaluate the effectiveness of closed incision negative pressure therapy (ciNPT) versus standard of care (SOC) dressings in reducing surgical site complications (SSCs).

Methods

- ► A total of 294 revision total knee arthroplasty (rTKA) patients (15 centres) at high risk for wound complications were randomised to ciNPT or SOC (n=147 each) and stratified by revision type (aseptic vs. septic). Demographics, comorbidities, causes of revision and duration of treatment were similar between cohorts (p>0.05).
- ► 242 patients with incisions completed follow-up, including 124 patients treated with Prevena Therapy (ciNPT) and 118 patients treated with an antimicrobial silver-impregnated dressing (SOC).
- ► Primary outcome was the 90-day incidence of SSCs with stratification in accordance with revision type. Secondary outcomes were the 90-day health care utilisation parameters (readmission, reoperation, dressing changes, and visits) and patient-reported outcomes (PRO). Treatment-related adverse events were compared and stratified as severe and non-severe.

Results

Compared to SOC, patients in the Prevena Therapy group demonstrated:

- Significantly decreased rates of surgical site complications (ciNPT 3.4" vs. SOC 14.3", p=0.0013*)
- Significantly lower readmission rates (ciNPT 3.4" vs. SOC 10.2", p=0.0208*)
- Reduced dressing changes (ciNPT 1.1±0.29 vs. SOC 1.3 ± 0.96 , p=0.0003*)

Conclusions

Prevena Therapy significantly mitigated 90-day surgical site complications, readmission rates, and reduced frequency of dressing changes compared with the standard of care among high-risk rTKA patients.

- Treatment-related adverse effects were similar between both cohorts.
- The benefit of ciNPT on specific SSCs and post-rTKA patient-reported outcomes (PRO) was not established and further studies are warranted.

(continued)

Calculation(s) are derived based on relative patient group incidence rate reported in this study.

*Statistically significant (p<0.05).

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com.

PRM in orthopaedic surgery PROMISES study Data

(continued)

Cost effectiveness

All patients:

\$989

Reduction in per-patient cost of care

\$1,047 3M™ Prevena™ Therapy vs. \$2,036 SOC

Higher-risk patients (CCI ≥2):

\$2,318

Reduction in per-patient cost of care

\$894 3M™ Prevena™ Therapy vs. \$3,212 SOC



Calculation(s) are derived based on relative patient group incidence rate reported in Cooper HJ, Bongards C, Silverman RP. Cost-Effectiveness of Closed Incision Negative Pressure Therapy for Surgical Site Management After Revision Total Knee Arthroplasty: Secondary Analysis of a Randomized Clinical Trial. *J Arthroplasty*. 2022 Aug;37(8S):S790–S795.

Read the full study here



Journal: The Journal of Arthroplasty

Title: The Effectiveness of Closed-Incision Negative-Pressure Therapy Versus Silver-Impregnated Dressings in Mitigating Surgical Site Complications in High-Risk Patients After Revision Knee Arthroplasty: The PROMISES Randomised Controlled Trial

Published: 5 March 2021

Higuera-Rueda CA, Emara AK, Nieves-Malloure Y, Klika AK, Cooper HJ, Cross MB, Guild GN, Nam D, Nett MP, Scuderi GR, Cushner FD, Piuzzi NS, Silverman RP. The Effectiveness of Closed-Incision Negative-Pressure Therapy Versus Silver-Impregnated Dressings in Mitigating Surgical Site Complications in High-Risk Patients After Revision Knee Arthroplasty: The PROMISES Randomised Controlled Trial. *J Arthroplasty*. 2021 Jul;36(7S):S295-S302.e14. **OPEN ACCESS** Note that the length of therapy may be outside the range recommended in the Instructions for Use.

^{*}Statistically significant (p<0.05).

PRM in orthopaedic surgery Anatone study

A risk-stratification algorithm to reduce superficial surgical site complications in primary hip and knee arthroplasty.

Anatone AJ, Shah RP, Jennings EL, Geller JA, Cooper HJ. A risk-stratification algorithm to reduce superficial surgical site complications in primary hip and knee arthroplasty. Arthroplasty Today. 2018 Dec;4(4):493-98.

Study design

Single institution retrospective review of records (United States).

Study purpose

The purpose of the Anatone study was to evaluate when to use 3M™ Prevena™ Therapy in primary total joint arthroplasties (TJAs). The author's risk stratification can be used as a potential guideline to identify patients that may benefit from Prevena Therapy.

Methods


- ▶ Patients were considered low risk if their calculated risk score was <2 and patients were considered high risk if their risk score was ≥2.
- A study population of 323 consecutive primary TJAs were evaluated, where 123 (38^u) of those patients were considered at elevated risk to receive Prevena Therapy. The remaining 200 patients received the standard postop dressing (AQUACEL® Ag SURGICAL cover dressing).
- ► A historical control population of 643 patients was identified who all received the standard postop dressing to test the impact of this risk score.
- ► Skin closure procedure was the same in both groups, and dressings were applied under sterile conditions in the operating room at the conclusion of the surgical procedure.
- ► The primary outcome measure was any postoperative surgical site complication (SSC†) that required intervention during the initial 90-day postoperative period.

Risk stratification algorithm scoring system

Risk factor	Weight	
BMI		
<18.5 kg/m²	1	
18.5-29.9 kg/m²	0	
30-34.9 kg/m ²	1	
35-39.9 kg/m ²	2	
>40 kg/m²	3	

Weight
2
1.3
1
1
2

Results



Guidance

The authors' risk stratification can be used as a potential guideline to identify patients who may benefit from Prevena Therapy.

Key points

- ► Among high-risk patients, there was a marked improvement in the rate of SSCs when treated prophylactically with Prevena Therapy as compared with historical controls (26.2" vs. 7.3"; p < 0.001).*
- ► Compared with historical controls, a modest but significant improvement in superficial SSCs after implementation of risk-stratification (12.0^u vs 6.8^u; p = 0.013) was observed.*
- ► Low-risk patients who continued to be treated with standard postop dressings in historical controls demonstrated no significant improvement $(8.6^{\rm u} \text{ vs } 6.5^{\rm u}; p = 0.344).$

Calculation(s) are derived based on relative patient group incidence rate reported in this study.

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com.

Read the full study here



Journal: Arthroplasty Today

Title: A risk-stratification algorithm to reduce superficial surgical site complications in primary hip and knee arthroplasty

Published: December 2018

Anatone AJ, Shah RP, Jennings EL, Geller JA, Cooper HJ. A risk-stratification algorithm to reduce superficial surgical site complications in primary hip and knee arthroplasty. Arthroplasty Today. 2018 Dec;4(4):493-98. OPEN ACCESS

^{*}Statistically significant (p<0.05).

[†]SSC was defined as any dehiscence, suture granuloma, drainage occurring beyond postoperative day 5, significant haematoma formation, or SSI as defined by the CDC that required unplanned postoperative interventions.

PRM in orthopaedic surgery Doman study

Comparison of surgical site complications with negative pressure wound therapy vs silver impregnated dressing in high-risk total knee arthroplasty patients: a matched cohort study.

Doman DM, Young AM, Buller LT, Deckard ER, Meneghini RM. Comparison of Surgical Site Complications With Negative Pressure Wound Therapy vs Silver Impregnated Dressing in High-Risk Total Knee Arthroplasty Patients: A Matched Cohort study. Journal of Arthroplasty. 2021 Oct;36(10):3437-3442.

Study design

Retrospective comparative cohort study (United States).


Study purpose

To compare high-risk primary TKA patients' rate of incisional and non-incisional wound complications, periprosthetic joint infections and reoperations.

Methods

- ► The 3M[™] Prevena[™] Therapy group comprised of 130 patients who had primary TKA between July 2018 and December 2019.
- ► The retrospective historical control group (AQUACEL® Ag SURGICAL) consisted of 130 patients, propensity matched 1:1, who underwent surgery between December 2016 and June 2018.
- High-risk criteria included active tobacco use, diabetes mellitus, BMI > 35 kg/m², autoimmune disease, chronic kidney disease, Staphylococcus aureus nasal colonisation, and non-aspirin anticoagulation.
- Study endpoints included incisional wound complications, defined as: cellulitis, focal swelling, suture reaction, dehiscence and Haematoma. Non-incisional wound complications were also assessed and defined as dressing reactions, blistering and rashes.

Calculation(s) are derived based on relative patient group incidence rate



Key points

- Among high-risk patients undergoing primary TKA, patients receiving Prevena Therapy had significantly fewer incisional wound complications when compared to patients receiving silver impregnated dressings.
- ► Although an increase in dressing reactions for Prevena Therapy patients was observed, the clinical impact was minimal.
- ► Results support the use of ciNPT as part of a risk mitigation strategy to reduce post operative complications in primary TKA.

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com.

Read the full study here



reported in this study.

Journal: The Journal of Arthroplasty

Title: Comparison of Surgical Site Complications With Negative Pressure Wound Therapy vs Silver Impregnated Dressing in High-Risk Total Knee Arthroplasty Patients: A Matched Cohort study

Published: 24 May 2021

Doman DM, Young AM, Buller LT, Deckard ER, Meneghini RM. Comparison of Surgical Site Complications With Negative Pressure Wound Therapy vs Silver Impregnated Dressing in High-Risk Total Knee Arthroplasty Patients: A Matched Cohort study. Journal of Arthroplasty. 2021 Oct;36(10):3437-3442. PMID 34140207.

^{*}Statistically significant (p<0.05).

PRM in orthopaedic surgery

Newman study

Use of closed incisional negative pressure wound therapy after revision total hip and knee arthroplasty in patients at high risk for infection: a prospective, randomised clinical trial.

Newman JM, Siqueira MBP, Klika AK, Molloy RM, Barsoum WK, Higuera CA. Use of Closed Incisional Negative Pressure Wound Therapy After Revision Total Hip and Knee Arthroplasty in Patients at High Risk for Infection: A Prospective, Randomised Clinical Trial. *Journal of Arthroplasty*. 2019 Mar;34(3):554–559.

Study design

Prospective, single-centre, randomised controlled trial (United States).

Study purpose


The purpose of the Newman study was to compare the use of Prevena Therapy to a sterile antimicrobial dressing (AQUACEL® Ag SURGICAL cover dressing) in revision arthroplasty (rTHA, rTKA) patients at high risk to develop wound complications.

Methods

- ► 160 patients undergoing elective rTHA and rTKA were prospectively randomised to receive Prevena Therapy or AQUACEL® Ag at a single institution.
- Patients had at least one risk factor for developing a wound complication.
- ► All patients received perioperative treatment and antibiotics.
- Study endpoints included wound complications (SSC including: SSIs, drainage and cellulitis), readmission and reoperation rates.
- ▶ Data collected at 2, 4 and 12 weeks postoperatively.

Key points

- High-risk patients could benefit from closed incision negative pressure therapy (ciNPT) to help reduce the risk of wound complications and reoperations after rTHA and rTKA.
- ► The authors suggest future multicentre clinical trials to further strengthen the results as well as a cost-benefit analysis.



(continued)

 $\label{lem:calculation} \mbox{Calculation(s) are derived based on relative patient group incidence rate reported in this study.}$

Although the authors reported use of Prevena Therapy for a mean of 3.6 days (ranging from 2 to 15 days), this mean time of application is outside the recommendations for Optimum Use as stated in the Prevena Incision Management System Clinician Guide Instructions for Use: The Prevena Incision Management System is to be continuously applied for a minimum of two days up to a maximum of seven days. Use for greater than 7 days is not recommended or promoted by 3M.

^{*}Statistically significant (p<0.05).

PRM in orthopaedic surgery Newman study

(continued)

Illustration of the 3M™ Prevena™ Therapy Incision Management System cost effectiveness based on Newman et al outcomes

Revision TKA surgery in high-risk population hypothetical economic model	Prevena Therapy	AQUACEL® Ag SURGICAL
Patients	79	80
Number of surgical site infections (a)	2	7
Cost per SSI ¹ (b)	€29,053	€29,053
Per patient infection cost (a*b)/n	€736	€2,542
Per patient therapy cost*	€295	€3
Total cost per patient	€1,031	€2,545
Potential per incision savings using Prevena Therapy	€	1.514

Potential cost savings Reduction in per patient cost for SSI €1,031 Prevena Therapy vs. €2,545 SOC

The infection cost assumption calculated form Hardstock et al. 2020 by subtracting the cost of a non-infected patient (13,781€) from the cost of an infected patient (42,834€) utilising 365-d follow-up costs (€) per patient-year.

*3M™ Prevena™ Peel and Place System Kit and AQUACEL® Ag SURGICAL price are an estimate; individual prices may vary.

The above model uses selected study data to provide an illustration of estimates of costs for use of the Prevena Therapy or AQUACEL® Ag SURGICAL. This model is an illustration and not a guarantee of actual individual costs, savings, outcomes or results. The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing.

Reference

1. Hardtstock, F., Heinrich, K., Wilke, T. et al. Burden of Staphylococcus aureus infections after orthopedic surgery in Germany. BMC Infect Dis 20, 233 (2020). https://doi.org/10.1186/s12879-020-04953-4 SSi cost calculated using 365-d follow-up costs (€) per patient-year for infected and non-infected patients.

Read the full study here



Journal: The Journal of Arthroplasty

Title: Use of Closed Incisional Negative Pressure Wound Therapy After Revision Total Hip and Knee Arthroplasty in Patients at High Risk for Infection: A Prospective, Randomised Clinical Trial

Published: 16 November 2018

Newman JM, Siqueira MBP, Klika AK, Molloy RM, Barsoum WK, Higuera CA. Use of Closed Incisional Negative Pressure Wound Therapy After Revision Total Hip and Knee Arthroplasty in Patients at High Risk for Infection: A Prospective, Randomised Clinical Trial. Journal of Arthroplasty. 2019 Mar;34(3):554-559. **OPEN ACCESS** Note that the length of therapy may be outside the range recommended in the Instructions for Use.

PRM in orthopaedic surgery Redfern study

Closed incision negative pressure therapy effects on postoperative infection and surgical site complication after total hip and knee arthroplasty.

Redfern RE, Cameron-Ruetz C, O'Drobinak SK, Chen JT, Beer KJ. Closed Incision Negative Pressure Therapy Effects on Postoperative Infection and Surgical Site Complication After Total Hip and Knee Arthroplasty. *J Arthroplasty.* 2017 Nov;32(11):3333–3339. doi: 10.1016/j.arth.2017.06.019. Epub 2017 Jun 17.

Study design

Single-centre, prospective versus historic control comparative study (United States).

Study purpose


The purpose of the Redfern study was to examine the use of closed incision negative pressure therapy (ciNPT) over clean closed surgical incisions after primary total joint replacement and whether 3M™ Prevena™ Therapy would reduce the rates of wound complications.

Methods

- ► The Prevena Therapy group was comprised of 192 patients representing 196 incisions, who were actively enrolled from 2013 to 2014.
- ► The historical control group consisted of 400 patients who underwent surgery from 2011 to 2012.
- Prevena Therapy was applied over the closed incision for 6-8 days postoperatively. The control group standard of care included a sterile gauze dressing with standard dressing changes.
- ► The rate of surgical site complications requiring medical or surgical intervention, including surgical site infections (deep and superficial infections), wound dehiscence, Haematomas, seromas, oedema/swelling, and drainage were compared between groups.

Key points

In this study, Prevena Therapy reduced the overall incidence of complications requiring medical or surgical intervention for hip and knee arthroplasty.



(continued)

Calculation(s) are derived based on relative patient group incidence rate reported in this study.

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at https://ncbgregulatory.3M.com.

^{*}Statistically significant (p<0.05).

PRM in orthopaedic surgery Redfern study

(continued)

Illustration of the 3M™ Prevena™ Therapy Incision Management System cost-effectiveness based on Redfern et al outcomes

Primary TKA/THA not limited to high-risk patients hypothetical economic model	Prevena Therapy	SOC – gauze dressing
Patients	196	400
Number of surgical site infections (a)	2	14
Cost per SSI¹ (b)	€29,053	€29,053
Per patient infection cost (a*b)/n	€296	€1,017
Per patient therapy cost*	€295	€3
Total cost per patient	€591	€1,020
Potential Per Incision Savings Using Prevena Therapy		€429



The infection cost assumption calculated form Hardstock et al. 2020 by subtracting the cost of a non-infected patient (13,781€) from the cost of an infected patient (42,834€) utilising 365-d follow-up costs (€) per patient-year.

*3M™ Prevena™ Peel and Place System Kit is an estimate; individual prices may vary.

The above model uses selected study data to provide an illustration of estimates of costs for use of the Prevena Therapy or gauze dressing. This model is an illustration and not a guarantee of actual individual costs, savings, outcomes or results. The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing.

Reference

1. Hardtstock, F., Heinrich, K., Wilke, T. et al. Burden of Staphylococcus aureus infections after orthopedic surgery in Germany. BMC Infect Dis 20, 233 (2020). https://doi.org/10.1186/s12879-020-04953-4 SSi cost calculated using 365-d follow-up costs (€) per patient-year for infected and non-infected patients.

Read the full study here



Journal: The Journal of Arthroplasty

Title: Closed Incision Negative Pressure Therapy Effects on Postoperative Infection and Surgical Site Complication After Total Hip and Knee Arthroplasty

Published: 16 June 2017

Redfern RE, Cameron-Ruetz C, O'Drobinak S, Chen J, Beer KJ. Closed incision negative pressure therapy effects on postoperative infection and surgical site complication after total hip and knee arthroplasty. Journal of Arthroplasty. 2017 Nov;32(11):3333-3339.

PMID 28705547 Note that the length of therapy may be outside the range recommended in the Instructions for Use.

PRM in orthopaedic surgery Cooper study

Randomised controlled trial of incisional negative pressure following high-risk direct total hip arthroplasty.

Cooper HJ, Santos WM, Neuwirth AL, Geller JA, Rodriguez JA, Rodriguez-Elizalde S, Shah RP. Randomised Controlled Trial of Incisional Negative Pressure Following High-Risk Direct Anterior Total Hip Arthroplasty. J Arthroplasty. 2022 Aug;37(8S):S931–S936.

Study type

This was a prospective randomised controlled trial (United States).

Study purpose

The purpose of this study is to determine whether ciNPT could decrease SSCs in high-risk patients undergoing DA THA. The direct anterior (DA) approach to total hip arthroplasty (THA) is associated with higher rates of surgical site complications (SSCs) compared to other approaches. Closed incision negative pressure therapy (ciNPT) is effective in reducing SSCs and surgical site infections (SSIs) in other populations.

Methods

- ► Population: study enrolled high-risk DA THA patients at 3 centres. Inclusion criteria was if subjects had previously identified risk factors for SSC: Body mass index (BMI) >30 kg/m², diabetes, active smoking or before hip surgery.
- ► Treatment: Patients were randomised after closure to either an occlusive (control) dressing or ciNPT dressing (3M™ Prevena™ Incision Management System) for 7 days. Both dressings were designed for 7 day use per manufacturer instructions.
- ► Follow up: All patients were followed for 90 days to assess SSCs.

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com.

Results

One hundred and twenty-two patients were enrolled and 120 completed the data collection. SSCs occurred in 18.3" (11/60) of control patients compared to 8.3^u (5/60) of ciNPT patients ($x^2 = 2.60$, P = .107).

- ► SSCs included dehiscence to the subcutaneous level (13) and prolonged drainage (3).
- ► Nine control (15.0^u) and 2 ciNPT (3.3^u) patients met CDC criteria for superficial SSI (P = .027).
- ► Fifteen of 16 SSCs resolved with local wound care. One in the ciNPT group required reoperation for acute PJI.

Conclusion

It was determined that among high-risk patients undergoing DA THA, there were lower rates of SSC and a significant reduction in the risk of superficial SSI with ciNPT.

Read the full study here



Journal: The Journal of Arthroplasty

Title: Randomised Controlled Trial of Incisional Negative Pressure Following High-Risk Direct Anterior Total Hip Arthroplasty

Published: 15 March 2022

Cooper HJ, Santos WM, Neuwirth AL, Geller JA, Rodriguez JA, Rodriguez-Elizalde S, Shah RP. Randomised Controlled Trial of Incisional Negative Pressure Following High-Risk Direct Anterior Total Hip Arthroplasty. J Arthroplasty. 2022 Aug;37(8S):S931-S936. doi: 10.1016/j. arth.2022.03.039. OPEN ACCESS

PRM in orthopaedic surgery Phillips study

Incisional negative pressure wound therapy in orthopaedic trauma: indications and outcomes.

Phillips, Rachel MD; Stannard, James P. MD; Crist, Brett D. MD. Incisional Negative Pressure Wound Therapy in Orthopaedic Trauma: Indications & Outcomes. Journal of Orthopaedic Trauma 36():p S22-S25, September 2022.

Study type

This was a literature review.

Study purpose

This review aims to discuss the indications and outcomes associated with the use of incisional negative pressure wound therapy (iNPWT) for the management of surgical incisions.

Outcomes

Indication for iNPWT: In patient population at high risk for developing SSIs, management of the surgical incision with iNPWT have reduced the incidence of SSIs.

Several meta-analyses and randomised controlled trials were evaluated to assess the efficacy of surgical site infections, wound dehiscence and other postoperative wound complications.

A 2019 meta-analysis analysed a total of 6 studies including 2 randomised controlled trials (RCTs) and 4 cohort studies comparing a mix of iNPWT systems to conventional wound dressings for closed incisions in orthopaedic trauma surgery found that 14 statistically significant lower incidence of deep SSIs (P = 0.002), superficial SSI (P = 0.03) and wound dehiscence (P = 0.02) was found in surgical incisions managed with iNPWT.

The results of 2 RCTs also support the use of iNPWT after primary and revision total joint arthroplasty. Total knee arthroplasty patients with a body mass index >35 kg/m² who were treated with incisional NPWT experienced fewer overall complications $(1.3^{\rm u} \text{ vs. } 21.6^{\rm u}; P = 0.01)$ and fewer dressing-related concerns (1.3 $^{\text{u}}$ vs. 10.8 $^{\text{u}}$; P = 0.01) compared with standard of care dressings.

Duration of treatment

Most studies that have reported the use of iNPWT before the availability of a portable device typically used iNPWT for 3–5 days during the inpatient hospital stay. More recent studies have extended therapy to 7 days. However, there are some contraindications to the iNPWT which includes if there is necrotic tissue with eschar present, preexisting infection, patients at high risk of excessive postoperative bleeding, and those who have an allergic reaction to any part of the NPWT system.

Conclusion

The literature review suggested that iNPWT seems to be an effective tool for decreasing the rates of surgical site infections and wound dehiscence across multiple specialties. SSI risk factors should be considered for either patients or wounds that are at high risk for infection and/or dehiscence.

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com.

Read the full study here



Journal: Journal of Orthopaedic Trauma

Title: Incisional Negative Pressure Wound Therapy in Orthopaedic Trauma: Indications & Outcomes

Published: September 2022

Phillips, Rachel MD; Stannard, James P. MD; Crist, Brett D. MD. Incisional Negative Pressure Wound Therapy in Orthopaedic Trauma: Indications & Outcomes. Journal of Orthopaedic Trauma 36():p S22-S25, September 2022. | DOI: 10.1097/ BOT.000000000002425 **OPEN ACCESS** PRM in orthopaedic surgery Zelle study

How Can negative pressure wound therapy pay for itself? – reducing complications is important.

Zelle BA, Kore L. How Can Negative Pressure Wound Therapy Pay for Itself? – Reducing Complications Is Important. *J Orthop Trauma*. 2022 Sep 1;36(Suppl 4):S31–S35.

Study type

This was a retrospective cohort study performed at a single, level-1 trauma centre using data from a lower extremity fracture registry (United States).

Study purpose

The purpose of this study was to investigate cost savings in high-risk fractures and to determine if the use of iNPWT (3M™ Prevena™ Therapy) in high-risk orthopaedic trauma patients reduces the costs. The hypothesis was that the use of iNPWT will provide an economic benefit in patients with OTA/AO type 41C and 43C closed fractures undergoing ORIF.

Methods

- Material: Patient data from single institution registry were retrospectively retrieved from January 2019 and September 2020.
- Population: The evaluation included all patients with closed OTA/AO type 41C or 43C fractures treated with ORIF (staged or immediately) during the study period.
- ▶ Procedure: Registry data were summarised to determine SSI rates in all patients with closed OTA/ AO type 41C and 43C fractures. 3 health economic models were developed using SSI rates of 13^u, 15^u and 17^u as reference rates. The incremental cost due to SSI was estimated to be \$51,364.

Result

Out of a total of 79 patients who underwent ORIF of a closed OTA/AO type 41C or 43C fractures, 27 (34^u) were deemed high risk for SSI and had iNPWT applied over the closed incision.

- ► There was no significant difference in rates of SSI when comparing iNPWT with non-iNPWT group (7.4" vs. 11.5", P = 0.7086).
- ► Patients in iNPWT group had the external fixator in place for a significantly longer time (10.6 days vs. 6.8 days; P = 0.0332). Length of hospital stay was longer for patients in the non-iNPWT group compared with the iNPWT group (10.2 vs. 5.4 days; P = 0.0155).
- ► Health economic models: For assumed SSI rates of 13^u, 15^u, and 17^u, the total infection costs for 100 patients would be \$667,732, \$770,460, and \$873,188, respectively, the per patient cost would be \$6,677, \$7,704, and \$8,732 respectively and iNPWT cohort, the total infection cost for 100 patients would be \$380,094 or \$3,801 per patient. Thus, when comparing the SSI rates, the differences in infection costs per patient were estimated to be \$2,381, \$3,409, and \$4,436, respectively. Hence, this health economic model suggests the use of the iNPWT in patients with high-risk OTA/AO type 41C and 43C fractures may provide estimated cost savings per patient that range between \$2,381 to \$4,436.

Conclusion

Based on this health economic model, the use of iNPWT (Prevena Therapy) may reduce the costs of SSI in high-risk orthopaedic trauma patients undergoing ORIF of their closed OTA/AO type 41C and 43C fractures.

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at https://ncbgregulatory.3M.com.

Read the full study here



Journal: Journal of Orthopaedic Trauma

Title: How Can Negative Pressure Wound Therapy Pay for Itself? – Reducing Complications Is Important

Published: September 2022


Zelle BA, Kore L. How Can Negative Pressure Wound Therapy Pay for Itself? – Reducing Complications Is Important. *J Orthop Trauma*. 2022 Sep 1;36(Suppl 4):S31–S35. doi: 10.1097/BOT.00000000000002427. **PMID: 35994307**.

PRM in orthopaedic surgery Decision guide

Decision guide

Patient and procedure risk stratification in orthopaedic surgery backed by clinical evidence.

While most surgical patients may benefit from 3M™ Prevena™ Therapy, patients at high risk for complications such as surgical site infection may see added benefit. The following uses select study data¹⁻³ to provide an illustrative guide to aid in risk stratification. This is not an all-inclusive list of risk factors. Clinicians are advised to use their clinical judgment to identify high-risk patients or high-risk procedures.



The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com.

1. Willy C, Agarwal A, Andersen CA, De Santis G, Gabriel A, Grauhan O, Guerra OM, Lipsky BA, Malas MB, Mathiesen LL, Singh DP, Reddy VS. Closed incision negative pressure therapy: international multidisciplinary consensus recommendations. Int Wound J. 2017 Apr;14(2):385-398. OPEN ACCESS 2. Higuera-Rueda CA, Emara AK, Nieves-Malloure Y, Klika AK, Cooper HJ, Cross MB, Guild GN, Nam D, Nett MP, Scuderi GR, Cushner FD, Piuzzi NS, Silverman RP. The Effectiveness of Closed-Incision Negative-Pressure Therapy Versus Silver-Impregnated Dressings in Mitigating Surgical Site Complications

in High-Risk Patients After Revision Knee Arthroplasty: The PROMISES Randomised Controlled Trial. J Arthroplasty. 2021 Jul;36(7S):S295-S302. e14. OPEN ACCESS Note that the length of therapy may be outside the range recommended in the Instructions for Use.

3. Newman JM, Siqueira MBP, Klika AK, Molloy RM, Barsoum WK, Higuera CA. Use of Closed Incisional Negative Pressure Wound Therapy After Revision Total Hip and Knee Arthroplasty in Patients at High Risk for Infection: A Prospective, Randomised Clinical Trial. Journal of Arthroplasty. 2019 Mar;34(3):554-559. OPEN ACCESS Note that the length of therapy may be outside the range recommended in the Instructions for Use.

Author biographies*



PRM in orthopaedic surgery Author biography - Cooper



H. John Cooper, MD

Associate Professor of Orthopaedic Surgery Columbia University Irving Medical Centre New York-Presbyterian Hospital, New York City, NY

Dr. Cooper is a paid consultant for 3M.

Originally from South Carolina, Dr. Cooper graduated from Duke University with a degree in mechanical engineering and materials science. He completed his medical education at Columbia University and his Orthopaedic residency at Lenox Hill Hospital, before spending a year in Chicago for a fellowship in adult reconstructive surgery at Rush University Medical Centre.

Dr. Cooper currently works as an associate professor of Orthopaedic surgery at Columbia University Irving Medical Centre in New York City. He has considerable experience in direct anterior hip arthroplasty, robotic knee arthroplasty, and complex primary and revision joint replacement.

Dr. Cooper is a well-respected clinician, educator and researcher. He has published over 130 peer-reviewed articles and book chapters on clinical outcomes and complications of hip and knee replacements and has been an invited and awarded speaker on these topics at national and international Orthopaedic meetings.

"I employ 3M™ Prevena™ Therapy as a proactive risk management tool, using an evidence-based approach to stratify patients on their unique patient-specific and procedure-specific risk factors. In my experience, proactively using Prevena Therapy on the high-risk patients has significantly improved their clinical outcomes (and mine as well)."

Dr. Cooper



Carlos Higuera-Rueda, MD

Cleveland Clinic Florida, Weston, FL

Dr. Higuera-Rueda is a paid consultant for 3M.

Dr. Carlos Higuera-Rueda is currently a staff surgeon at the Cleveland Clinic Florida, where he divides his time between leadership, research and patient care. He is the Chairman of the Levitetz Department of Orthopaedic Surgery at Cleveland Clinic Florida and Director of the Orthopaedic and Rheumatology Centre. Dr. Higuera completed his residency at the Cleveland Clinic and a clinical fellowship at Thomas Jefferson University Hospital.

Dr. Higuera specialises in hip and knee arthroplasty surgery. He uses alternative approaches for primary hip and knee arthroplasty to optimise recovery. He is interested in complex revision procedures including infections. His research interest is mainly in periprosthetic joint infections including diagnostic tools, patient optimissation and overall outcomes after arthroplasty. He is currently working on developing new technologies to diagnose and treat such infections. He is the past-president of the Musculoskeletal Infection Society.

"Based on the level 1 clinical evidence in adult reconstruction revision surgery, we use 3M™ Prevena™ Therapy on our high-risk patient population to reduce the risk of SSC, SSI, readmissions and reoperations. In our experience, the portability and ease-of-use of the technology has also helped to reduce length of stay and office visits."

Dr. Higuera-Rueda

PRM in orthopaedic surgery Author biography - Crist



Brett D. Crist, MD, **FACS, FAAOS**

Professor Vice Chair of Business Development Director Orthopaedic Trauma Service Director Orthopaedic Trauma Fellowship Department of Orthopaedic Surgery University of Missouri School of Medicine, Columbia, MO

Dr. Crist is a paid consultant for 3M.

After obtaining a bachelor's degree from Tabor College in Hillsboro, Kansas, Dr. Crist earned his medical degree from the University of Kansas School of Medicine. He completed his residency at the University of Kansas School of Medicine, Wichita, and a fellowship in Orthopaedic trauma at the University of California-Davis.

Dr. Crist specialises in Orthopaedic trauma/fracture care, limb deformity correction, hip and pelvis reconstruction including total hip arthroplasty, and young adult hip disorders/hip preservation. Areas of interest include:

- Anterior total hip arthroplasty
- ▶ Fractures
- ► Hip and pelvic reconstruction Surgery
- ► Hip arthroscopy
- ► Minimally invasive surgery

- ► Orthopaedic rehabilitation
- Orthopaedic trauma surgery
- ► Pelvic surgery
- ► Skeletal trauma
- ► Limb deformity correction

"In my practice, I have standardised my approach for using 3M™ Prevena™ Therapy. Leveraging Proactive Risk Management (PRM), I stratify my patients based on common procedural/ patient risk factors to reduce the risk of SSIs, thereby improving patient outcomes. I place a Prevena dressing on most of my high-risk patients."

Dr. Crist

PRM in plastic surgery

3M.de/Prevena 3M™ Prevena Restor™ Dressings can be used on a variety of anatomical locations.

Patients that used 3M™ Prevena™ Therapy experienced reduced complications and reoperation after breast reconstruction.

Gabriel A, et al. (Loma Linda University). The Impact of Closed Incision Negative Pressure Therapy on Postoperative Breast Reconstruction Outcomes. *Plast Reconstr Surg Glob Open.* 2018;6:e1880.

Study design

Retrospective, comparative study (United States).

Study purpose

The investigators compared incision management outcomes in patients who received Prevena Therapy versus standard of care (SOC) after breast reconstruction mastectomy.

Methods

- ► Single site retrospective observational study: 2009–2017.
- ▶ 356 patients (Prevena Therapy n=177 v SOC n=179).
- ► 665 closed breast incisions (Prevena Therapy n=331 vs. SOC n=334).
- ► SOC: 3M[™] Steri-Strip[™] Wound Closures.
- ► 3M[™] Prevena[™] Plus Customizable Dressing.
- ► Patients were discharged home after 1 night stay and returned for follow-up on POD 3 and 7.
- ► Patient demographics, chemotherapy exposure, surgical technique, number of drains, time to drain removal, and 90-day postoperative complication rates were analysed.

Summary of findings

The use of Prevena Therapy following post-mastectomy breast reconstruction was associated with significantly lower rates of infection, dehiscence, necrosis and seromas. A significantly shorter time to drain removal and fewer returns to the OR were also achieved.

In addition to the above observed clinical outcomes, an economic analysis relying on this study data showed a mean per patient cost saving for SSC of \$218.1

\$2,010 Prevena Therapy vs. \$2,228 standard of care.

Results Reduction in SSCs* Reduction in SSIs* 8.5^u (28/331) Prevena Therapy vs. 2.1^u (7/331) Prevena Therapy vs. 4.5" (15/334) SOC (p=0.0225)* 15.9" (53/334) SOC (p=0.0092)* Reduction in seroma* Reduction in reoperations* 1.8" (6/331) Prevena Therapy vs. 2.4" (8/331) Prevena Therapy vs. 5.7" (19/334) SOC (p=0.0106)* 5.4" (18/334) SOC (p=0.0496)* Reduction in dehiscence* Reduction in necrosis*1 2.4^u (8/331) Prevena Therapy vs. 5.1^u (17/331) Prevena Therapy vs. 5.4" (18/334) SOC (p=0.0178)*1 9.3" (31/334) SOC (p=0.0070)*1

(continued)

Cost assessment includes variable hospital costs (for both the index hospitalisation and all readmission days within 30 days related to any wound complication). Hospital variable costs (not charges) for each admission were obtained from hospital administration.

Calculation(s) are derived based on relative patient group incidence rate reported in this study. *Statistically significant (p<0.05).

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at https://ncbgregulatory.3M.com.

References

1. Gabriel A, Maxwell P. Economic analysis based on the use of closed-incision negative-pressure therapy after postoperative breast reconstruction. Plast Reconstr Surg 2019;143:36S. PRM in plastic surgery Ferrando study

(continued)

Read the full study here



Journal: Plastic and Reconstructive Surgery – Global Open

Title: The Impact of Closed Incision Negative Pressure Therapy on Postoperative Breast **Reconstruction Outcomes**

Published: August 2018

Gabriel A, Sigalove S, Sigalove N, Storm-Dickerson T, Rice J, Maxwell P, Griffin L. The impact of closed incision negative pressure therapy on postoperative breast reconstruction outcomes. Plastic and Reconstructive Surgery Global Open. 2018 Aug; 6(8):e1880. **OPEN ACCESS**

PRM in plastic surgery Ferrando study

Improved outcomes with the use of 3M™ Prevena™ Therapy after breast surgery in high risk patients.

Ferrando PM, Ala A, Bussone R, Bergamasco L, Actis Perinetti F, Malan F. Closed Incision Negative Pressure Therapy in Oncological Breast Surgery: Comparison with Standard Care Dressings. *Plast Reconstr Surg Glob Open.* 2018 Jun 15;6(6):e1732.

Study design

Prospective, comparative study (Italy).

Study purpose

Evaluated the use of Prevena Therapy for oncological breast surgery patients that were high-risk for unfavourable healing.


Methods

- ► From January 2015 to June 2015, 47 patients were prospectively selected. Patients were undergoing oncological breast surgery.
- ► Inclusion criteria: patients had a minimum of 4 risk factors with at least 1 high risk factor.
- ▶ 17 patients (25 surgeries) voluntary treated with ciNPT; the remaining 20 patients (22 surgeries) chose conventional post-surgery dressing.
- ► SOC: 3M[™] Steri-Strip[™] Wound Closures.
- ► 3M[™] Prevena[™] Plus Customizable Dressing for 7 days.
- ▶ 90 days follow-up to evaluate postsurgical complications.
- ► At 12 months, the quality of life, scar, and overall aesthetic outcomes were assessed.

Summary of study findings

This study demonstrates that the use of Prevena Therapy in oncological breast surgery resulted in a statistically significant reduction in surgical site complications.

At the 12-month follow-up, questionnaires completed by both the plastic surgeon (Observer Scar Assessment Scale) and the patient (Patient Scar Assessment Scale) on level of satisfaction showed a significant difference in favour of Prevena Therapy.



(continued)

Calculation(s) are derived based on relative patient group incidence rate reported in this study.


^{*}Statistically significant (p<0.05).

PRM in plastic surgery Savage study

(continued)

Illustration of the 3M™ Prevena™ Therapy Incision Management System cost-effectiveness based on Ferrando et al outcomes

Oncological breast surgery hypothetical economic model	Prevena Therapy	Steri-strip
Patients	25	22
Number of surgical site infections (a)	1	10
Cost per SSI¹ (b)	€6,133	€6,133
Per patient infection cost (a*b)/n	€245	€2,788
Per patient therapy cost*	€385	-
Total cost per patient	€630	€2,788
Potential per incision savings using Prevena Therapy		€2,158



Hypothetical attributable cost of GER SSI/SSC taken from DRG Code J25Z Cancer related simple skin sparing mastectomy for low risk patient, no NPWT used, no post operative complications €3,464.63: Compared to J06Z. Cancer related complex breast reconstruction pathway for high risk patient, use of NPWT and 14 days LOS €9,598.09. Resulting in hypothetic cost of SSI/SSC as €6,133 (Complex Mamma Wound Care Grouping J35Z only 0.01% or 147 cases per year reached diagnosis DRG: J35Z 21 Days 14,400 €).

*3M™ Prevena™ Peel and Place System Kit is an estimate; individual prices may vary.

The above model uses selected study data to provide an illustration of estimates of costs for use of the Prevena Therapy or gauze dressing. This model is an illustration and not a guarantee of actual individual costs, savings, outcomes or results. The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing.

Reference

1. https://www.g-drg.de/content/download/10834/file/ Fallpauschalenkatalog_2022_20211123.pdf

Read the full study here



Journal: Plastic and Reconstruction Surgery - Global Open

Title: Closed Incision Negative Pressure Therapy in Oncological Breast Surgery: Comparison with Standard Care Dressings

Published: 15 June 2018

Ferrando PM, Ala A, Bussone R, Bergamasco L, Actis Perinetti F, Malan F. Closed Incision Negative Pressure Therapy in Oncological Breast Surgery: Comparison with Standard Care Dressings. Plast Reconstr Surg Glob Open. 2018 Jun 15;6(6):e1732. OPEN ACCESS PRM in plastic surgery Savage study

Reduced wound complications and opioid use with the use of 3M™ Prevena™ Therapy after bilateral breast reduction.

Savage N, Jain M, Champion R, Snell B. Incisional negative pressure wound therapy in bilateral breast reduction patients. *Australas J Plast Surg.* 2020; 3(1):30–38.

Study design

Retrospective comparative cohort study (Australia).

Study purpose


Evaluate the effect of Prevena Therapy, on surgical complications, opioid use and hospitalisation length after bilateral breast reduction.

Methods

- ► Consecutive bilateral breast reductions performed by a single surgeon June 2015 to August 2017. 52 patients analyzed: SOC (n=29) and Prevena Therapy (n=23).
- ► Prevena Therapy was used for 7 days with no drains and no fitted garment.
- ► SOC: application of an adhesive non-woven fabric dressing, gauze and adhesive fabric dressing again, drains removed on post-operative day 1, fitted garment used post OP.
- Discharge criteria defined as able to mobilise, subjective pain score less than 4, feeling subjectively well.
- Outcome Measure: SSC including local inflammatory response, dehiscence, surgical site infection, delayed healing, nipple necrosis, abscess; opioid use measured in oral morphine equivalents.

Summary of study findings

- ► This is the first study to provide evidence for the use of ciNPT in bilateral breast reduction. This study indicates that Prevena Therapy could be associated with a significant reduction in surgical site complication occurrences, decreased total ward opioid use and discharge opioid prescription as well as decreased hospital length of stay.
- ► The study was not limited to high-risk patients.



(continued)

Calculation(s) are derived based on relative patient group incidence rate reported in this study.


^{*}Statistically significant (p<0.05).

PRM in plastic surgery Ayuso study

(continued)

Illustration of the 3M™ Prevena™ Therapy Incision Management System cost-effectiveness based on Savage et al outcomes

Oncological breast surgery hypothetical economic model	Prevena Therapy	Steri-strip
Patients	23	29
Number of surgical site infections (a)	3	13
Cost per SSI¹ (b)	€6,133	€6,133
Per patient infection cost (a*b)/n	€800	€2,749
Per patient therapy cost*	€385	=
Total cost per patient	€1,185	€2,749
Potential per incision savings using Prevena Therapy		€1,564



Hypothetical attributable cost of GER SSI/SSC taken from DRG Code J25Z Cancer related simple skin sparing mastectomy for low risk patient, no NPWT $\,$ used, no post operative complications €3,464.63: Compared to J06Z. Cancer related complex breast reconstruction pathway for high risk patient, use of NPWT and 14 days LOS \leqslant 9,598.09. Resulting in hypothetic cost of SSI/SSC as €6,133 (Complex Mamma Wound Care Grouping J35Z only 0.01% or 147 cases per year reached diagnosis DRG: J35Z 21 Days 14,400 €).

*3M™ Prevena™ Peel and Place System Kit is an estimate; individual prices may vary.

The above model uses selected study data to provide an illustration of estimates of costs for use of the Prevena Therapy or gauze dressing. This model is an $\,$ illustration and not a guarantee of actual individual costs, savings, outcomes or results. The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing.

Reference

1. https://www.g-drg.de/content/download/10834/file/ Fallpauschalenkatalog_2022_20211123.pdf

Read the full study here



Journal: Australian Journal of Plastic Surgery

Title: Incisional negative pressure wound therapy in bilateral breast reduction patients

Published: 23 March 2020

Savage N, Jain M, Champion R, Snell B. Incisional negative pressure wound therapy in bilateral breast reduction patients. Australas J Plast Surg. 2020; 3(1):30-38.

PRM in plastic surgery Ayuso study

ciNPT for open abdominal wall reconstruction with concomitant panniculectomy.

Ayuso SA, Elhage SA, Okorji LM, Kercher KW, Colavita PD, Heniford BT, Augenstein VA. Closed-Incision Negative Pressure Therapy Decreases Wound Morbidity in Open Abdominal Wall Reconstruction With Concomitant Panniculectomy. *Ann Plast Surg.* 2022 Apr 1;88(4):429–433.

Study design

Retrospective cohort study (United States).

Study purpose

To evaluate the use of closed-incision negative pressure therapy (ciNPT) and its effects on postoperative wound complications in open Abdominal Wall Reconstruction (AWR) patients with Concomitant Panniculectomy (CP).

Methods


- Prospective institutional database identified 67 patients that received 3M™ Prevena™ Therapy. These patients were matched 1:1 to 67 patients that received standard surgical dressings before the use of ciNPT.
- ► In the study period, patient prehabilitation and perioperative protocols at the institution were the same which aids in eliminating confounders.
- ► From 2016 onward all patient rehabilitation and perioperative protocols at the institution were the same
- Prevena Therapy was used for 7 days.
- Concomitant Panniculectomy makes this a study on high-risk patients.
- Primary outcomes: wound complications defined as seroma requiring drainage, cellulitis requiring antibiotics, deep wound infection and superficial wound breakdown.

Key points

Patients undergoing abdominal wall reconstruction with concomitant panniculectomy can be at higher risk for wound complications due to the need for large incisions and tissue undermining. In this study, the use of Prevena Therapy helped significantly decrease the risk of postoperative wound occurrences including superficial wound breakdown. The study also demonstrated the lessened need for wound-related reoperations in ciNPT patients.

Calculation(s) are derived based on relative patient group incidence rate reported in this study.

*Statistically significant (p<0.05)




(continued)

PRM in plastic surgery Savage study

(continued)

Illustration of the 3M™ Prevena™ Therapy Incision Management System cost-effectiveness based on Ayuso et al outcomes

Plastics AWR with CP hypothetical economic model	Prevena Therapy	Control
Number of patients (n)	67	67
Number of surgical site infections (a)	10	24
Cost per SSI¹ (b)	€11,545	€11,545
Per patient infection cost (a*b)/n	€1,723	€4,315
Per patient therapy cost*	€295	_
Total cost per patient	€2,018	€4,315
Potential per incision savings using Prevena Therapy		€2,297



The above model uses selected study data to provide an illustration of estimates of costs for use of the Prevena Therapy or gauze dressing. This model is an illustration and not a guarantee of actual individual costs, savings, outcomes or results. The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing

*3M™ Prevena™ Peel and Place System Kit is an estimate; individual prices may vary.

Attributable cost of SSI in abdominal surgery obtained from the median difference in resource costs including inpatient costs, surgery costs, surgical ward costs, medication costs, laboratory and diagnostic costs of €11,545 illustrated in Strobel et al. 2020.

Reference

1. Strobel, R.M., Leonhardt, M., Förster, F. et al. The impact of surgical site infection—a cost analysis. *Langenbecks Arch Surg* 407, 819–828 (2022). https://doi.org/10.1007/s00423-021-02346-y.

Read the full study here



Journal: Annals of Plastic Surgery

Title: Closed-Incision Negative Pressure Therapy Decreases Wound Morbidity in Open Abdominal Wall Reconstruction With Concomitant Panniculectomy

Published: April 2022

Ayuso SA, Elhage SA, Okorji LM, Kercher KW, Colavita PD, Heniford BT, Augenstein VA. Closed-Incision Negative Pressure Therapy Decreases Wound Morbidity in Open Abdominal Wall Reconstruction With Concomitant Panniculectomy. Ann Plast Surg. 2022 Apr 1;88(4):429-433. doi: 10.1097/ SAP.0000000000002966. **PMID: 34670966.**

PRM in plastic surgery Pieri study

Managing no-drain mastectomy with closed incision negative pressure wound therapy using full-coverage foam dressings.

Pieri A, Aisling E, Kay K, Irving J, Robert T, Cain H, Kalra L, Critchley A. Managing No-Drain Mastectomy with Closed Incision Negative Pressure Wound Therapy using Full-Coverage Foam Dressings. *European Journal of Surgical Oncology.* 2023 Feb; 49(2):e94.

Study Design

This was a single centre, case-control trial (United Kingdom).

Study Purpose

The trial was aimed to evaluate whether the mastectomies managed with ciNPT using full coverage foam dressings exhibited reduced need for seroma intervention and reduced seroma aspiration volumes.

Methods

Seroma intervention data was retrospectively gathered from a single centre for patients undergoing simple mastectomy, mastectomy with sentinel lymph node biopsy, or mastectomy with axillary lymph node clearance. 30 sequential patients treated with conventional dressings in control arm and 25 sequential patients treated with ciNPT with full-coverage foam dressings (3M™ Prevena Restor™ Bella·Form™ Dressing) were selected for intervention arm.

Results

There were 31 mastectomy cases in each arm (including bilateral cases). There was no significant difference in surgery type between the groups.

- 1. Compared to control group, fewer patients in the intervention group developed postoperative seroma (20 control versus 15 intervention).
- 2. More subjects needed aspiration in control group than intervention group (16 control vs 12 intervention).
- 3. Fewer visits to the seroma clinic were needed for intervention group than control group (1 control vs. 0 intervention, p=0.012).
- 4. Intervention group had lower total aspiration volumes (843ml control vs. 368ml intervention, p=0.023).

Conclusion

The study indicated that the patients managed with ciNPT with full-coverage foam dressings required fewer seroma-related clinical episodes and experienced reduced total seroma volume. The use of ciNPT has reduced the costs and improved the services and therefore it has been adopted as the standard practice at this centre.

Read the full study here



Journal: European Journal of Surgical Oncology

Title: Managing No-Drain Mastectomy with Closed Incision Negative Pressure Wound Therapy using Full-Coverage Foam Dressings

Published: February 2023

Pieri A, Aisling E, Kay K, Irving J, Robert T, Cain H, Kalra L, Critchley A. Managing No-Drain Mastectomy with Closed Incision Negative Pressure Wound Therapy using Full-Coverage Foam Dressings. *European Journal of Surgical Oncology.* 2023 Feb; 49(2):e94. doi: 10.1016/j.ejso.2022.11.287.

PRM in plastic surgery Silverman study

The use of closed incision negative pressure therapy for incision and surrounding soft tissue management: Expert panel consensus recommendations.

Ronald P. Silverman MD, John Apostolides MD, FACS, Abhishek Chatterjee MD, MBA, Anthony N. Dardano DO, FACS, Regina M. Fearmonti MD, FACS Allen Gabriel MD, FACS Robert T. Grant MD, MSc, FACS, Owen N. Johnson III MD, FACS, Suresh Koneru MD, Anna A. Kuang MD, Andrea A. Moreira MD, Steven R. Sigalove MD, FACS

Study type

The study type was an Expert Panel convened to develop consensus recommendations. In the absence of high-quality studies, an expert panel of plastic surgeons reviewed the current literature and formed consensus utilising a modified Delphi technique.

Study purpose

The purpose of the study was to identify conditions in which ciNPT with full-coverage dressings is most appropriate, and address challenges to the implementation and sustainability of ciNPT.

Methods

Consensus building was done using modified Delphi technique, which involved three rounds of input to gather feedback and identify topics with potential for agreement. Consensus was defined as ≥80^u agreement among panel members.

Selected panelists had experience using ciNPT with both conventional and novel dressings, previously presented or published on the use of ciNPT, were able to present their cases demonstrating use of ciNPT in the panel meetings and were able to understand and participate in consensus formation process.

The panel recommended use of ciNPT with full-coverage dressings when 2 or more risk factors for surgical site complications are present.

Results

The panel was able to establish 10 consensus statements. Recommendations for the use of ciNPT with full coverage dressings were provided for patient and incision related risk factors, therapy duration, appropriate pressure settings to be used, and lastly, techniques used for ciNPT. The panel recommended that future studies on ciNPT should focus on identifying the benefits of use and overcoming implementation barriers.

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com.

Read the full study here



Journal: International Wound Journal

Title: The use of closed incision negative pressure therapy for incision and surrounding soft tissue management: Expert panel consensus recommendations

Published: 21 August 2021

Silverman, RP, Apostolides, J, Chatterjee, A, et al. The use of closed incision negative pressure therapy for incision and surrounding soft tissue management: Expert panel consensus recommendations. Int Wound J. 2022;19(3): 643-655. https://doi.org/10.1111/iwj.13662


OPEN ACCESS

PRM in plastic surgery Decision guide

Decision guide

Patient and procedure risk stratification in plastic surgery backed by clinical evidence.

While most surgical patients may benefit from 3M™ Prevena™ Therapy, patients at high risk for complications such as surgical site infection may see added benefit. The following uses select study data¹⁻² to provide an illustrative guide to aid in risk stratification. This is not an all-inclusive list of risk factors. Clinicians are advised to use their clinical judgment to identify high-risk patients or high-risk procedures.



The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com.

1. Willy C, Agarwal A, Andersen CA, De Santis G, Gabriel A, Grauhan O, Guerra OM, Lipsky BA, Malas MB, Mathiesen LL, Singh DP, Reddy VS. Closed incision negative pressure therapy: international multidisciplinary consensus recommendations. Int Wound J. 2017 Apr;14(2):385-398. OPEN ACCESS 2. Gabriel A, Sigalove S, Sigalove N, Storm-Dickerson T, Rice J, Maxwell P, Griffin L. The impact of closed incision negative pressure therapy on postoperative breast reconstruction outcomes. Plastic and Reconstructive Surgery Global Open. 2018 Aug; 6(8):e1880. OPEN ACCESS

Author biographies*



PRM in plastic surgery Author biography - Chatterjee



Abhishek Chatterjee, MD, MBA

Chief of the Division of Plastic Surgery Division of Surgical Oncology Tufts Medical Centre Boston, MA

Dr. Chatterjee is a paid consultant for 3M.

Dr. Chatterjee is a board-certified plastic surgery and fellowship trained breast oncologic surgeon practicing at Tufts Medical Centre in Boston, MA. After completing his MD/MBA training at the University of Connecticut, he went on to do eight years of surgical residency at Dartmouth Hitchcock Medical Centre in New Hampshire in plastic surgery and followed this with a one-year breast surgical oncology fellowship at the University of Pennsylvania. With this unique training in oncology and plastic surgery, much of Dr. Chatterjee's practice involves the removal of cancer and the reconstruction using oncoplastic surgical techniques. He is active within his own institution as the President of the Medical Staff and sits on several committees as a member in both national breast oncologic and plastic surgery societies. He is presently Associate Professor of Surgery at Tufts Medical Centre and is the Chief of Plastic Surgery.

Academically, he enjoys training surgical residents daily and has published more than 90 peer-reviewed journal articles, most of which are either first or senior authored.

"My use of ciNPT began when I wanted to reduce my wound complication rates in high-risk breast cancer patients, so that I could get my patients to adjuvant therapy after surgery without delay. Now I continue to use ciNPT on all of my patients with any high-risk incisions to decrease my overall complication rates regardless of anatomical location."

Dr. Chatterjee

PRM in plastic surgery Author biography - Gabriel



Allen Gabriel, MD, **FACS**

Private Practice Vancouver, Washington

Dr. Gabriel is a paid consultant for 3M.

Allen Gabriel, MD, is an Assistant Professor and Director of Research in the Department of Surgery at Loma Linda University, Loma Linda, California. He is a board-certified plastic surgeon that believes plastic and reconstructive surgery provides a unique opportunity to deal with a wide variety of needs ranging from addressing congenital anomalies, to breast reconstruction following mastectomy, to aesthetic procedures such as breast and facial cosmetic procedures.

In 2001, Dr. Gabriel was chosen by the prestigious Loma Linda University to join the Integrated Plastic Surgery Residency Program. While at Loma Linda University, he volunteered on a medical mission to Ethiopia with Operation Good Samaritan. In addition, he served on several leadership committees and was the chief resident prior to completing his residency. In 2007, Dr. Gabriel was selected by Dr. G. Patrick Maxwell to enter a Breast and Aesthetic Surgery Fellowship in conjunction with Baptist Hospital in Nashville, Tennessee. Completion of this program provided him with advanced training in breast and aesthetic surgery.

Dr. Gabriel is one of the few medical students in the country to have received the prestigious Humanism in Medicine Award. This award led to the creation of the University of Nevada's Humanism in Medicine Honor Society, of which Dr. Gabriel is still an active member. During medical school, he was involved with both clinical and basic science research, earning several research awards and publications prior to graduating. Dr. Gabriel has been invited to speak nationally and internationally on breast and aesthetic surgery. Dr. Gabriel is a Fellow of the American College of Surgeons. He is also a member of several prestigious organisations including the American Board of Plastic Surgery, American Society of Bariatric Plastic Surgeons, American Society of Plastic Surgeons, and California Society of Plastic Surgeons.

Since 1995, Dr. Gabriel has authored more than three dozen abstracts and chapters in peer-reviewed publications, including articles on liposuction, tummy tuck, breast anatomy and breast embryology.

"In 2012, we started using closed incision negative pressure therapy in complex reconstructions in my practice. Subsequently in 2014, we decided to expand use of the technology into breast reconstructions because of the positive clinical results on key patient outcomes. At that time, my colleagues wanted to better understand how to leverage a risk stratification algorithm to inform a more standardised approach of the therapy. We then published the figure, which we still use today."

Dr. Gabriel


(continued)

PRM in plastic surgery

Author biography – Gabriel

(continued)

Incisions at risk for surgical complications¹



BMI – body mass index; DM – diabetes mellitus

Checklist of potential risk factors for surgical complications.

Reference

1. Gabriel A, Sigalove SR, Maxwell GP. Initial Experience Using Closed Incision Negative Pressure Therapy after Immediate Postmastectomy Breast Reconstruction. *Plastic and Reconstructive Surgery Global Open.* 2016 Jul 22;4(7):e819. **OPEN ACCESS**

PRM in plastic surgery

Author biography – Pieri



Andrew Pieri MBBS, MRes, FRCS

The Newcastle upon Tyne Hospitals NHS Foundation Trust · Department of Surgery

Mr Pieri is a paid consultant for 3M.

Mr Pieri is a consultant oncoplastic breast surgeon working within the NHS at the Royal Victoria Infirmary in Newcastle upon Tyne and at several private hospitals in the Northeast of England, performing cancer and cosmetic breast surgery. He performs a comprehensive range of oncoplastic procedures including perforator flaps, mammoplasty surgery and autologous or implant-based breast reconstructions. In the private sector, he performs both cancer and cosmetic breast surgery.

With numerous research publications and a surgical device patent to his name, Mr Pieri has a keen interest in surgical device innovation. He has introduced and formally evaluated a number of novel technologies. His unit was the first in the UK to progress from wire-guided breast conservation surgery – Mr Pieri published the benefits of seeds over wires. He has recently introduced ICG-guided sentinel node biopsy and axilla reverse mapping and is working with industry to develop a Newcastle-based programme of education for colleagues in this technique.

He is education lead for the breast unit in Newcastle, being involved in both undergraduate and post graduate education programmes.

"The question is no longer whether you should use negative pressure therapy to reduce the risk of surgical site complications; it's in which patients and in which procedures should it be used.

We have been using Prevena Therapy for many years in patients undergoing therapeutic mammoplasties, implant reconstructions or patients with high risk factors such as obesity or active smoking."

Mr. Pieri



Patients that used 3M™ Prevena™ Therapy experienced significant reduction in complications, reoperation and readmission rates for high-risk groin incision procedures.

Kwon J, Staley C, McCullough M, Goss S, Arosemena M, Abai B, Salvatore D, Reiter D, DiMuzio P. A randomised clinical trial evaluating negative pressure therapy to decrease vascular groin incision complications. J Vasc Surg. 2018 Dec;68(6):1744-1752.

Study design

Prospective, single-centre, randomised controlled trial (United States).

Study purpose

This prospective RCT evaluated negative pressure therapy (Prevena Therapy) to decrease wound complications and associated healthcare costs.

Methods

- ► The study included 119 femoral incisions closed primarily after elective vascular surgery procedures.
- High-risk inclusion criteria: re-operative surgery, BMI > 30, pannus, prosthetic graft, poor nutrition, immunosuppression, or HbA1c>8.
- ▶ 1:1 Randomised to standard gauze (n=60) vs. Prevena Therapy (n=59).
- Outcomes evaluated at post-operative day 30: SSI, wound complications, length of stay (LOS), reoperation, readmission.

Summary of findings

Study suggests that negative pressure therapy for patients at high risk for groin wound complications:

- ► Significantly reduces major wound complication
- ► Significantly reduces reoperation and readmission rates
- ► Closed incision negative pressure therapy (ciNPT) may lead to a reduction in hospital cost

Calculation(s) are derived based on relative patient group incidence rate reported in this study. *Statistically significant (p<0.05).

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com

Results Reduction in SSCs* Reduction in SSIs* 11.9^u (7/59) Prevena Therapy vs. 10.1^u (6/59) Prevena Therapy vs. 26.7" (16/60) SOC (p=0.001)* 31.6" (12/60) SOC (p=0.001) Reduction in return to OR* Reduction in readmissions* 8.5^u (5/59) Prevena Therapy vs. 6.8^u (4/59) Prevena Therapy vs. 18.3" (11/60) SOC (p=0.05)* 16.7" (10/60) SOC (p=0.04)*

ciNPT is recommended for all groin incisions considered at high risk for wound complications. In addition to the above observed clinical outcomes, this study data¹ showed per patient cost saving of \$6,045 for Prevena Therapy patients.

\$30,492 Prevena Therapy vs. \$36,537 SOC

Cost assessment includes variable hospital costs (for both the index hospitalisation and all readmission days within 30 days related to any wound complication). Hospital variable costs (not charges) for each admission were obtained from hospital administration.

Reference

1. Cost Assessment includes variable hospital costs (for both the index hospitalisation and all readmission days within 30 days related to any wound complication). Hospital variable costs (not charges) for each admission were obtained from hospital administration.

Read the full study here



Journal: Journal of Vascular Surgery Title: A randomised clinical trial evaluating negative pressure therapy to decrease vascular groin incision complications

Published: 17 August 2018

Kwon J, Staley C, McCullough M, Goss S, Arosemena M, Abai B, Salvatore D, Reiter D, DiMuzio P. A randomised clinical trial evaluating negative pressure therapy to decrease vascular groin incision complications. J Vasc Surg. 2018 Dec;68(6): 1744-1752. doi: 10.1016/j.jvs.2018.05.224.

OPEN ACCESS

PRM in vascular surgery Pleger study

Reduction of groin wound complications in vascular surgery patients using closed incision negative pressure therapy (ciNPT): a prospective, randomised, single-institution study.

Pleger SP, Nink N, Elzien M, Kunold A, Koshty A, Böning A. Reduction of groin wound complications in vascular surgery patients using closed incision negative pressure therapy (ciNPT): a prospective, randomised, single-institution study. *International Wound Journal*. 2018 Feb;15(1):75–83. **OPEN ACCESS**

Study design

Single centre randomised controlled trial (Germany).

Study purpose


The purpose of the study was to investigate the effectiveness of ciNPT (3M™ Prevena™ Therapy) compared to conventional therapy on vascular surgical groin incisions.

Methods

- Patients were randomised and treated with either Prevena Therapy or the control therapy, a conventional adhesive plaster.
- ► 100 patients with 129 groin incisions were analysed: ciNPT consisted of 58 incisions; Control consisted of 71 incisions.
- ► Inclusion criteria for high-risk patients: age > 50 years, diabetes mellitus, renal insufficiency, malnutrition, obesity and chronic obstructive pulmonary disease.
- ciNPT was applied intraoperatively and removed on days 5–7 postoperatively.
- ► Wound evaluation based on the Szilagyi classification took place postoperatively on days 5–7 and 30.

Key points

This study found that the use of ciNPT demonstrated a statistically significant reduction of postoperative wound healing complications in the groin on postoperative days 5–7 and 30-day revision surgery.



(continued)

 $\label{eq:Calculation} \mbox{(s) are derived based on relative patient group incidence rate reported in this study}$

Although the authors reported use of Prevena Therapy for a mean of 3.6 days (ranging from 2 to 15 days), this mean time of application is outside the recommendations for Optimum Use as stated in the 3M™ Prevena™ Incision Management System Clinician Guide Instructions for Use: The Prevena Incision Management System is to be continuously applied for a minimum of two days up to a maximum of seven days. Use for greater than 7 days is not recommended or promoted by 3M.


^{*}Statistically significant (p<0.05)

PRM in vascular surgery Pleger study

(continued)

Illustration of the 3M™ Prevena™ Therapy Incision Management System cost-effectiveness based on Pleger et al outcomes

Vascular groin hypothetical economic model	Prevena Therapy	Control
Number of incisions (n)	58	71
Number of surgical site infections (a)	5	30
Cost per infection ¹ (b)	€10,262	€10,262
Cost of infection per incision (a*b)/n	€884	€4,336
Cost of therapy per incision*	€295	-
Total cost per incision	€1,179	€4,336
Potential per incision savings using Prevena Therapy		€3,157



Hypothetical cost of SSI taken from DRG Code F08C Non –Diabetic Bypass Patient with SSI, Dehiscence, Revision Surgery €17,308.52 : Compared to F59C. Non-Diabetic Vascular DRG with mean LOS Post-Op Stay of 6.3 Days €7,046.52. Resulting in hypothetic cost of Ssi/Complications as **10,262.33**€

The above model uses selected study data to provide an illustration of estimates of costs for use of the Prevena Therapy or gauze dressing. This model is an illustration and not a guarantee of actual individual costs, savings, outcomes or results. The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing.

Reference

1. https://www.g-drg.de/content/download/10834/file/ Fallpauschalenkatalog_2022_20211123.pdf

Read the full study here



Journal: International Wound Journal

Title: Reduction of groin wound complications in vascular surgery patients using closed incision negative pressure therapy (ciNPT): a prospective, randomised, single-institution study

Published: 25 October 2017

Pleger SP, Nink N, Elzien M, Kunold A, Koshty A, Böning A. Reduction of groin wound complications in vascular surgery patients using closed incision negative pressure therapy (ciNPT): a prospective, randomised, single-institution study. International Wound Journal. 2018 Feb;15(1):75-83. **OPEN ACCESS**

PRM in vascular surgery Gombert study

Closed incision negative pressure therapy reduces surgical site infections in vascular surgery: a prospective randomised trial (AIMS trial).

Gombert A, Babilon M, Barbati ME, Keszei A, von Trotha KT, Jalaie H, Kalder J, Kotelis D, Greiner A, Langer S, Jacobs MJ, Grommes J. Closed Incision Negative Pressure Therapy Reduces Surgical Site Infections in Vascular Surgery: A Prospective Randomised Trial (AIMS Trial). *Eur J Vasc Endovasc Surg.* 2018 Sept; 56(3):442–448.

Study design

Prospective, multi-centre, randomised controlled trial (Germany).

Study purpose


This prospective RCT aimed to assess the potential benefit of ciNPT (3M™ Prevena™ Therapy) application to reduce the surgical site infection risk after groin incision for vascular surgery.

Methods

- ► The study evaluated 188 patients who underwent vascular surgery for peripheral artery disease (PAD) with a longitudinal groin incision at two sites in Germany between July 2015 and May 2017.
- ► High-risk inclusion criteria: smoking, cardiac risk factors including hypertension, coronary heart disease, or history of myocardial infarction, metabolic disorders including diabetes, dyslipidaemia, hyperhomocysteinaemia or chronic renal failure.
- When a groin incision was performed on both sides, only one side was randomised and assessed for the study.
- 30-day SSIs were assessed using the Szilagyi classification.

Key points

- Study found closed incision negative pressure therapy (ciNPT) was associated with a reduced incidence of SSIs when compared to control group.
- ► High-risk patients could benefit from ciNPT to help reduce the risk of total SSI.



(continued)

Calculation(s) are derived based on relative patient group incidence rate reported in this study.

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at https://ncbgregulatory.3M.com.


^{*}Statistically significant (p<0.05).

PRM in vascular surgery Gombert study

(continued)

Illustration of the 3M™ Prevena™ Therapy Incision Management System cost-effectiveness based on Gombert et al outcomes

Vascular groin hypothetical economic model	Prevena Therapy	Control
Number of patients (n)	98	90
Number of Surgical Site Infections (a)	13	30
Cost per SSI¹ (b)	€10,292	€10,292
Cost of SSI per patient (a*b)/n	€1,409	€3,420
Cost of therapy per patient*	€295	
Total cost per patient	€1,704	€3,420
Potential per incision savings using Prevena Therapy		€1,716



Hypothetical cost of SSI taken from DRG Code F08C Non -Diabetic Bypass Patient with SSI, Dehiscence, Revision Surgery € 17,308.52: Compared to F59C. Non-Diabetic Vascular DRG with mean LOS Post-Op Stay of 6.3 Days €7,046.52. Resulting in hypothetic cost of Ssi/Complications as €10,262.33

*3M™ Prevena™ Peel and Place System Kit is an estimate; individual prices may vary.

The above model uses selected study data to provide an illustration of estimates of costs for use of the Prevena Therapy. This model is an illustration and not a guarantee of actual individual costs, savings, outcomes or results. The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing.

1. https://www.g-drg.de/content/download/10834/file/ Fallpauschalenkatalog_2022_20211123.pdf

Read the full study here



Journal: European Journal of Vascular and Endovascular Surgery

Title: Closed Incision Negative Pressure Therapy Reduces Surgical Site Infections in Vascular Surgery: A Prospective Randomised Trial (AIMS Trial)

Published: 2 July 2018

Gombert A, Babilon M, Barbati ME, Keszei A, von Trotha KT, Jalaie H, Kalder J, Kotelis D, Greiner A, Langer S, Jacobs MJ, Grommes J. Closed Incision Negative Pressure Therapy Reduces Surgical Site Infections in Vascular Surgery: A Prospective Randomised Trial (AIMS Trial). Eur J Vasc Endovasc Surg. 2018 Sept; 56(3):442-448. OPEN ACCESS

PRM in vascular surgery Gombert, Dillavou study

A systematic review and meta-analysis of randomised controlled trials for the reduction of surgical site infection in closed incision management versus standard of care dressings over closed vascular groin incisions.

Gombert A, Dillavou E, D'Agostino R Jr, Griffin L, Robertson JM, Eells M. A systematic review and meta-analysis of randomised controlled trials for the reduction of surgical site infection in closed incision management versus standard of care dressings over closed vascular groin incisions. *Vascular*. 2020 Jun;28(3):274–284.

Study type

This study was a systematic review and meta-analysis of randomised controlled trials.

Study purpose

The purpose of the study was to assess the effect of ciNPT (3M[™] Prevena[™] Incision Management System; KCI, San Antonio, TX) versus traditional postsurgical dressing use on SSI rates over closed groin incisions following vascular surgery.

Methods

- ▶ Literature search: A systematic literature search using PubMed, OVID, EMBASE and QUOSA was performed on 3 January 2019 by two independent reviewers to assess the literature between 1 January 2005 and 31 December 2018. The review conformed to the statement and reporting checklist of the Preferred Reporting Items for Systematic Reviews and Meta-Analyses.
- ► Inclusion criteria: Inclusion criteria were abstracts or manuscripts written in English, published study, conference abstract, randomised controlled trial (RCT), comparison of ciNPT use over closed groin incisions to traditional postoperative dressings, endpoint/outcome of SSI, and a study population ≥10.
- ► Data type: Characteristics of study participants, surgical procedure, type of dressing used, duration of treatment, incidence of surgical site infection, and length of follow-up were extracted.
- ► Statistical methods used: The odds ratios (OR) were calculated to assess the effect of ciNPT versus SOC on vascular groin incision SSIs. Weighted odds ratios and 95^u confidence intervals (CI) were calculated to pool study and control groups in each publication for analysis. high-risk patients, normal-risk patients, and Szilagyi I, II, III outcomes were assessed between ciNPT and control groups.

Results

Out of 615 publications that were identified during the literature search, 303 abstracts and titles were screened against the inclusion and exclusion criteria. Six RCTs were included in the analysis. The screening process is shown in Figure 1. There was a total of 733 incisions, of which 362 (49.4^u) received ciNPT and 371 (50.6^u) received standard of care. Patients treated with ciNPT had a lower risk of developing an SSI when compared to the control arm (OR = 3.06, $95^{\rm u}$ CI [2.05, 4.58], p < 0.05) showing highly significant effect in favour of ciNPT. High-risk, normal-risk, Szilagyi I, and Szilagyi II meta-analyses were also statistically significant in favour of ciNPT use (p < 0.05). However, risk of bias in selecting meta-analysis, differences in inclusion/exclusion criteria and selection of procedure type pose major limitation for the study.




Figure 1: A Prisma flowchart showing the process of identifying articles to be included in systemic review and meta-analysis.

(continued)

PRM in vascular surgery Gombert, Dillavou study

(continued)

Conclusion

The study shows that ciNPT usage demonstrated a statistically significant reduction in the incidence of SSI relative to traditional postsurgical dressings in patients undergoing vascular surgery with groin incisions.

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com.

Read the full study here



Journal: Vascular

Title: A systematic review and meta-analysis of randomised controlled trials for the reduction of surgical site infection in closed incision management versus standard of care dressings over closed vascular groin incisions

Published: 19 January 2020

Gombert A, Dillavou E, D'Agostino R Jr, Griffin L, Robertson JM, Eells M. A systematic review and meta-analysis of randomised controlled trials for the reduction of surgical site infection in closed incision management versus standard of care dressings over closed vascular groin incisions. Vascular. 2020 Jun;28(3):274-284.

OPEN ACCESS

PRM in vascular surgery Bloom study

A cost-utility analysis of the use of closed-incision negative pressure system in vascular surgery groin incisions.

Bloom JA, Tian T, Homsy C, Singhal D, Salehi P, Chatterjee A. A Cost-Utility Analysis of the Use of Closed-Incision Negative Pressure System in Vascular Surgery Groin Incisions. The American Surgeon. 2022;0(0).

Study type

The study was literature review looking at prospective randomised control trials that determined the probabilities and outcomes for femoral-popliteal bypass with and without ciNPT.

Study purpose

The aim of the study was to perform a cost-effectiveness analysis evaluating closed incision negative pressure therapy (ciNPT, 3M™ Prevena™ Incision Management System, KCI Medical San Antonio, TX) use in femoral-popliteal bypass with prosthetic graft.

Methods

Population selected: 65-year-old male with Vascular surgery such as lower extremity claudication and tissue loss.

Model: The model used femoral-popliteal graft with vs without prosthetic graft. Under each decision tree data was obtained incorporating the probability of health states and the costs and utilities associated with them such as post operative minor and major wound infections, sartorius flap reconstruction, excision of graft and axillary femoral bypass, amputation and death.

Analysis

Data from retrospective analysis was used to create a Decision analysis tree to highlight the more

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in the all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com.

cost-effective strategy. Cost data from Medicare charge and reimbursement defined as sum of hospital cost and surgeon reimbursement fees were utilised. Utility scores converted to QALYs obtained for all health states from previously published utility scores representing health states ranging from 0 (death) to 1 (healthy) were used for analysis. Incremental cost effectiveness ratio (ICER) was performed with willingness to pay \$50,000.

Results

The decision tree analysis demonstrated that femoral-popliteal bypass with Prevena Therapy has a higher clinical effectiveness (QALY) of 6.14 compared to without Prevena Therapy (6.13) and is more cost effective (with \$40,138 vs without \$41,774) resulting in a negative ICER of -234,764.03, favouring ciNPT. This indicated a dominant strategy.

In one-way sensitivity analysis, femoral-popliteal bypass without Prevena Therapy was cost-effective strategy if the probability of successful surgery in the Prevena Therapy arm was less than 84.9^u or if cost of Prevena Therapy exceeds \$3,139.

Conclusion

Despite the added cost of Prevena Therapy, its use is more cost-effective in vascular surgical operations using groin incisions.

Read the full study here



Journal: The American Surgeon

Title: A Cost-Utility Analysis of the Use of Closed-Incision Negative Pressure System in Vascular Surgery Groin Incisions

Published: 7 April 2022

Bloom JA, Tian T, Homsy C, Singhal D, Salehi P, Chatterjee A. A Cost-Utility Analysis of the Use of Closed-Incision Negative Pressure System in Vascular Surgery Groin Incisions. The American Surgeon. 2022;0(0).


doi:10.1177/00031348221087395

Decision guide PRM in vascular surgery

Decision guide

Patient and procedure risk stratification in vascular surgery backed by clinical evidence.


While most surgical patients may benefit from 3M™ Prevena™ Therapy, patients at high risk for complications such as surgical site infection may see added benefit. The following uses select study data¹⁻⁴ to provide an illustrative guide to aid in risk stratification. This is not an all-inclusive list of risk factors. Clinicians are advised to use their clinical judgment to identify high-risk patients or high-risk procedures.



The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com.

1. Willy C, Agarwal A, Andersen CA, De Santis G, Gabriel A, Grauhan O, Guerra OM, Lipsky BA, Malas MB, Mathiesen LL, Singh DP, Reddy VS. Closed incision negative pressure therapy: international multidisciplinary consensus recommendations. Int Wound J. 2017 Apr;14(2):385-398. OPEN ACCESS 2. Kwon J, Staley C, McCullough M, Goss S, Arosemena M, Abai B, Salvatore D, Reiter D, DiMuzio P. A randomised clinical trial evaluating negative pressure therapy to decrease vascular groin incision complications. J Vasc Surg. 2018 Dec;68(6):1744-1752. OPEN ACCESS 3. Gombert A, Babilon M, Barbati ME, Keszei A, von Trotha KT, Jalaie H, Kalder J, Kotelis D, Greiner A, Langer S, Jacobs MJ, Grommes J. Closed Incision Negative Pressure Therapy Reduces Surgical Site Infections in Vascular Surgery: A Prospective Randomised Trial (AIMS Trial). Eur J Vasc Endovasc Surg. 2018 Sept; 56(3):442-448. OPEN ACCESS 4. Pleger SP, Nink N, Elzien M, Kunold A, Koshty A, Böning A. Reduction of groin wound complications in vascular surgery patients using closed incision negative pressure therapy (ciNPT): a prospective, randomised, single-institution study. International Wound Journal. 2018 Feb;15(1):75-83. **OPEN ACCESS**

Author biographies*



PRM in vascular surgery Author biography - Dillavou



Ellen Dillavou, MD, **FACS, RPVI**

Medical Director, Vascular Surgery WakeMed Hospitals Raleigh, NC

Dr. Dillavou is a paid consultant for 3M.

Ellen Dillavou, MD, FACS, is the medical director of vascular surgery at the WakeMed hospital system in Raleigh, NC. She earned a BA at Macalester College in St. Paul, MN, an MD at the University of Arizona, completed general surgery training at Thomas Jefferson University of Philadelphia, and a vascular surgery fellowship at The University of Pittsburgh Medical Centre. Her work centres on complicated dialysis access, surgical quality improvement and surgical site infection prevention.

"I became aware of 3M™ Prevena™ Therapy while investigating interventions that help mitigate the risk of Surgical Site Infections. As I dug into the research, it became guite clear that Prevena is one of the most impactful therapies available to reduce SSIs for groin incisions in vascular surgery. I now use Prevena on all of my patients who are considered high risk for incisions at the groin or below."

Dr. Dillavou

PRM in vascular surgery Author biography - Gombert



PD Dr. med. Alexander Gombert, PhD, FEBVS


Endovascular Specialist, Consultant of vascular surgery, European Centre of Vascular Surgery Aachen-Maastricht, Clinic for Vascular Surgery, University of Aachen, Germany

Dr. Gombert is a paid consultant for 3M.

PD Dr. med. Gombert was born in 1983. He is working as a consultant of Vascular Surgery at one of the largest centres for Vascular Surgery in Germany, the European Vascular Centre Aachen-Maastricht. He is the initiator and principal investigator of the "Aachen Incision management system (AIMS) trial," a randomised, prospective, multicentre study, comparing the effect of 3M™ Prevena™ incision management system with standard wound dressings after groin incision for vascular surgical procedures. Furthermore, he is establishing one of the biggest databases for tissue samples of patients undergoing thoracoabdominal aortic surgery. Beneath his activity in the fields of wound healing and thoracic aortic aneurysm research, he is working in the venous research group of the European Vascular Centre Aachen-Maastricht. He is an active reviewer of different high-ranked vascular surgery journals. PD Dr. Gombert is the author of several high-ranked peer-reviewed publications focusing on different aspects of vascular surgery. Furthermore, he is frequently invited to speak at vascular surgical and general surgical meetings around the world. He is living together with his wife and three children in the area of Aachen.

"Prevena Therapy is an extremely valuable proactive risk management tool that can help improve patient outcomes, while reducing costs associated with surgical site infections (SSIs). With more than 200 peer-reviewed publications studying Prevena, several common patient and procedural risk factors within the literature have been elevated to help support clinical decision making. In my practice, we utilise Prevena on every at-risk patient and procedure, advancing the standard of care for surgical patients."

Dr. Gombert



High risk, obese sternotomy patients that used 3M[™] Prevena[™] Therapy experienced significant reduction in rate of wound infection.

Grauhan O, Navasardyan A, Hofmann M et al. Prevention of post sternotomy wound infections in obese patients by negative pressure wound therapy. J Thorac Cardiovasc Surg 2013;145:1387–1392.

Study design

Prospective, single-centre, controlled trial (Germany).

Study purpose


To evaluate negative pressure wound dressing treatment (Prevena Therapy) for infection prevention.

Methods

- ► The study included 150 consecutive obese patients who underwent a median sternotomy at a single site in Germany between April 2010 and October 2011.
- ► Inclusion criteria was a body mass index ≥ 30kg/m².
- ► The control group, (conventional wound dressings) consisted of 75 patients. Post Op dressing change day 1–2.
- ciNPT (Prevena Therapy) group consisted of 75 patients. Placed immediately after suturing. Post Op dressing removal at day 6–7.
- ► The primary end point was wound infection within 90 days.

Summary of study findings

Closed incision negative pressure therapy (ciNPT) reduces the rate of post sternotomy wound infection in high-risk, obese patients.



(continued)

 $\label{lem:calculation} \mbox{Calculation(s) are derived based on relative patient group incidence rate reported in this study.}$


The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at https://ncbgregulatory.3M.com.

^{*}Statistically significant (p<0.05).

(continued)

Illustration of the 3M™ Prevena™ Therapy Incision Management System cost-effectiveness based on Grauhan et al (2013) outcomes

Hypothetical economic model	Prevena Therapy	Control
Number of patients (n)	75	75
Number of Surgical Site Infections (a)	3	12
Cost per SSI¹ (b)	€22,905	€22,905
Cost of SSI per patient (a*b)/n	€916	€3,664
Cost of therapy per patient*	€295	-
Total cost per patient	€1,211	€3,664
Potential per incision savings using Prevena Therapy		€2,453



Difference in the median cost of infected CABG patient (${\leqslant}36,261$) and non-infected pateint (\in 13,355) is an attributable median cost of SSI of \in 22,905.

*3M™ Prevena™ Peel and Place System Kit is an estimate; individual prices

The above model uses selected study data to provide an illustration of estimates of costs for use of the Prevena Therapy. This model is an illustration and not a guarantee of actual individual costs, savings, outcomes or results. The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing.

Reference

1. Graf K, Ott E, Vonberg RP, Kuehn C, Haverich A, Chaberny IF. Economic aspects of deep sternal wound infections. Eur J Cardiothorac Surg. 2010 Apr;37(4):893-6. doi: 10.1016/j.ejcts.2009.10.005. Epub 2009 Nov 6. PMID: 19896860.

Read the full study here



Journal: The Journal of Thoracic and Cardiovascular Surgery

Title: Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy

Published: 29 October 2012

Grauhan O, Navasardyan A, Hofmann M, Muller P, Stein J, Hetzer R. Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy. Journal of Thoracic and Cardiovascular Surgery. 2013 May;145(5):1387–92. **OPEN ACCESS**

Effect of surgical incision management on wound infections in post sternotomy patient population.

Grauhan O, Navasardyan A, Tutkun B *et al.* Effect of surgical incision management on wound infections in a post sternotomy patient population. *Int Wound J* 2014;11:6–9.

Study design

Prospective, single-centre, controlled trial (Germany).

Study purpose


To evaluate Prevena Therapy vs. conventional wound dressings over closed surgical incisions in reducing wound infections.

Methods

- ► The study group (Prevena Therapy) included <u>all</u> prospective patients undergoing median sternotomy from September–October 2013 totaling 237 patients.
- The control group (conventional wound dressings) included <u>all</u> median sternotomy patients retrospectively analysed for the period of January 2008 December 2009 totalling 3,508 patients.
- ► No defined High Risk Inclusion Criteria.
- Prevena Therapy placed immediately after suturing.
 Post Op dressing removal at day 6–7.
- ► The primary end point was wound infection within 30 days.

Summary of study findings

Application of surgical incision management using ciNPT on clean, closed surgical incisions reduced the rate of post sternotomy wound infection.



(continued)

Calculation(s) are derived based on relative patient group incidence rate reported in this study.


The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at https://ncbgregulatory.3M.com.

^{*}Statistically significant (p<0.05).

(continued)

Illustration of the 3M™ Prevena™ Therapy Incision Management System cost-effectiveness based on Grauhan et al (2014) outcomes

Hypothetical economic model	Prevena Therapy	Control
Number of patients (n)	237	3,508
Number of Surgical Site Infections (a)	3	119
Cost per SSI ¹ (b)	€22,905	€22,905
Cost of SSI per patient (a*b)/n	€289	€776
Cost of therapy per patient*	€295	-
Total cost per patient	€585	€776
Potential per incision savings using Prevena Therapy		€191



Difference in the median cost of infected CABG patient (€36,261) and non-infected pateint (\in 13,355) is an attributable median cost of SSI of **€22.905.**

*3M™ Prevena™ Peel and Place System Kit is an estimate; individual prices may vary.

The above model uses selected study data to provide an illustration of estimates of costs for use of the Prevena Therapy. This model is an illustration and not a guarantee of actual individual costs, savings, outcomes or results. The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing.

Reference

1. Graf K, Ott E, Vonberg RP, Kuehn C, Haverich A, Chaberny IF. Economic aspects of deep sternal wound infections. Eur J Cardiothorac Surg. 2010 Apr;37(4):893-6. doi:10.1016/j.ejcts.2009.10.005. Epub 2009 Nov 6. PMID: 19896860

Read the full study here



Journal: International Wound Journal

Title: Effect of surgical incision management on wound infections in a poststernotomy patient population

Published: 23 May 2014

Grauhan O, Navasardyan A, Tutkun B, Hennig F, Müller P, Hummel M, Hetzer R. Effect of surgical incision management on wound infections in a poststernotomy patient population. Int Wound J.

2014 Jun;11 Suppl 1(Suppl 1):6-9.

OPEN ACCESS

The impact of closed incision negative pressure therapy on prevention of median sternotomy infection for high risk cases: a single centre retrospective study.

Suelo-Calanao RL, Thomson R, Read M, Matheson E, Isaac E, Chaudhry M, Loubani M. The impact of closed incision negative pressure therapy on prevention of median sternotomy infection for high risk cases: a single centre retrospective study. *Journal of Cardiothoracic Surgery.* 2020 Aug 19;15(1):222. **OPEN ACCESS**

Study design

Retrospective cohort study (United Kingdom).

Study purpose


To assess the effect of closed incision negative pressure therapy (ciNPT) on the infection rate of patients at high risk for sternal wound infection (SWI).

Methods

- ► This study included patients who underwent full median sternotomies between January 2009 to December 2016.
- ► Retrospective study included patients 3 years before the introduction of ciNPT (3M[™] Prevena[™] Therapy) and 3 years after introduction.
- ► No clinician change in practice other than the use of Prevena Therapy for high-risk patients.
- High-Risk patients: ≥ 2 risk factors: Obesity, COPD, Age ≥ 80, Diabetes.
- All patients followed up at 6 weeks following discharge.
- ► Before the introduction of ciNPT, 162 high-risk patients received SOC. After the introduction of ciNPT, 158 received ciNPT.

Key point

ciNPT reduced the incidence of post sternotomy sternal wound infections (SWIs) in high-risk patients.



(continued)

 $\label{eq:calculation} \mbox{Calculation(s) are derived based on relative patient group incidence rate reported in this study.}$

^{*}Statistically significant (p<0.05).

PRM in cardiothoracic surgery Suelo-Calanao study

(continued)

Illustration of the 3M™ Prevena™ Therapy Incision Management System cost-effectiveness based on Suelo-Calanao et al 2020 clinical outcomes

Potential per incision savings using Prevena Therapy		€1,228
Total cost per patient	€1,599	€2,827
Cost of therapy Per patient*	€295	-
Cost of SSI per patient (a*b)/n	€1,305	€2,827
Cost per Surgical Site Infection ¹ (b)	€22,905	€22,905
Number of Surgical Site Infections (a)	9	20
Number of patients (n)	158	162
Hypothetical economic model	Prevena Therapy	Control



Difference in the median cost of infected CABG patient (€36,261) and non-infected pateint (€13,355) is an attributable median cost of SSI of €22,905.

*3M™ Prevena™ Peel and Place System Kit is an estimates; individual prices may vary.

The above model uses selected study data to provide an illustration of estimates of costs for use of the Prevena™ Therapy or Standard of Care (Control). This model is an illustration and not a quarantee of actual individual costs, savings, outcomes or results. Results are based on selected study data and may not be typical. The hospital is advised to use this model as an illustration only to assist in an overall assessment of products and pricing.

Reference

1. Graf K, Ott E, Vonberg RP, Kuehn C, Haverich A, Chaberny IF. Economic aspects of deep sternal wound infections. Eur J Cardiothorac Surg. 2010 Apr;37(4):893-6. doi: 10.1016/j.ejcts.2009.10.005. Epub 2009 Nov 6. PMID: 19896860.

Read the full study here



Journal: Journal of Cardiothoracic Surgery

Title: The impact of closed incision negative pressure therapy on prevention of median sternotomy infection for high risk cases: a single centre retrospective study


Published: 19 August 2020

Suelo-Calanao RL, Thomson R, Read M, Matheson E, Isaac E, Chaudhry M, Loubani M. The impact of closed incision negative pressure therapy on prevention of median sternotomy infection for high risk cases: a single centre retrospective study. Journal of Cardiothoracic Surgery. 2020 Aug 19;15(1):222. OPEN ACCESS PRM in sardiothoracic surgery Decision guide

Decision guide

Patient and procedure risk stratification in cardiothoracic surgery backed by clinical evidence.

While most surgical patients may benefit from 3M™ Prevena™ Therapy, patients at high risk for complications such as surgical site infection may see added benefit. The following uses select study data¹⁻³ to provide an illustrative guide to aid in risk stratification. This is not an all-inclusive list of risk factors. Clinicians are advised to use their clinical judgment to identify high-risk patients or high-risk procedures.



The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com.

1. Willy C, Agarwal A, Andersen CA, De Santis G, Gabriel A, Grauhan O, Guerra OM, Lipsky BA, Malas MB, Mathiesen LL, Singh DP, Reddy VS. Closed incision negative pressure therapy: international multidisciplinary consensus recommendations. Int Wound J. 2017 Apr;14(2):385-398. OPEN ACCESS

2. Grauhan O, Navasardyan A, Hofmann M, Muller P, Stein J, Hetzer R. Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy. Journal of Thoracic and Cardiovascular Surgery. 2013 May;145(5):1387-92. **OPEN ACCESS**

3. Suelo-Calanao RL, Thomson R, Read M, Matheson E, Isaac E, Chaudhry M, Loubani M. The impact of closed incision negative pressure therapy on prevention of median sternotomy infection for high risk cases: a single centre retrospective study. Journal of Cardiothoracic Surgery. 2020 Aug 19;15(1):222. **OPEN ACCESS**



Read the full Willy Consensus study here

Author biographies*





V. Sreenath (Seenu) Reddy, MD, MBA, FACS

Chief, Division of Cardiothoracic Surgery, Centennial Heart & Vascular, Nashville, TN

Dr. Reddy is a paid consultant for 3M.

V. Sreenath (Seenu) Reddy, MD, MBA, FACS is Chief, Division of Cardiothoracic Surgery at Centennial Heart & Vascular Centre in Nashville, TN. He earned his Medical Doctorate from The University of Alabama School of Medicine. He then served his internship and completed a residency in General Surgery at Vanderbilt University Medical Centre. Dr. Reddy then received his training in Cardiovascular and Thoracic Surgery at Emory University Medical Centre. In addition, he completed a fellowship in advanced endovascular surgery at Emory University Medical Centre.

"The available clinical evidence in vascular, plastic, orthopaedic, cardiothoracic and spine surgery demonstrates that 3M™ Prevena™ Therapy should be the standard of care for high-risk patients or high-risk procedures. We have integrated Proactive Risk Management, or PRM, into my practice and routinely use Prevena on these groups of patients."

Dr. Reddy





Negative pressure wound therapy for surgical-site infections: a randomised trial.

Javed A, Teinor J, Wright M, Ding D, Burkhart R, Hundt J, Cameron J, Makary M, He J, Eckhauser F, Wolfgang C, Weiss M. Negative pressure wound therapy for surgical-site infections: a randomised trial. *Annals of Surgery.* 2019; 269(6):1034–1040.

Study design

Randomised controlled trial, single-centre (John Hokins Hospital, United States).

Study purpose

The purpose of the Javed RCT was to evaluate efficacy of closed incision negative pressure therapy (ciNPT), Prevena Therapy, to decrease surgical site infections (SSI) after open pancreaticonduodenectomy.

Methods


- Patients undergoing pancreaticoduodenectomy procedures were eligible if considered to be high risk for SSI.
- ► High risk for SSI was defined as a risk score of ≥1 defined by Poruk et al[†] where preoperative bile stent/ drain received 1 point and neoadjuvant chemotherapy received 1 point. Points were summed for each patient.
- ► A total of 123 patients analysed: Prevena Therapy (n=62) v. Standard of Care (SOC) (n=61).
- Preoperative and operative characteristics were not significantly different between the two groups.
- ► The primary outcome was 30-day SSI (superficial or deep).

Summary of findings

This randomised controlled trial from Johns Hopkins Hospital demonstrated significantly lower SSI rates in high-risk patients receiving 3M™ Prevena™ Therapy after pancreaticoduodenectomy (31.1^u vs. 9.7^u; p=0.003)*.

SSIs resulted in an increased hospitalisation cost of **\$9,778 per patient.**

Implementing Prevena Therapy into surgical practice can help reduce the risk of potential complications and associated healthcare costs.



Calculation(s) are derived based on relative patient group incidence rate reported in this study.

- *Statistically significant (p<0.05).
- [†]Poruk et al. A novel, validated risk score to predict surgical site infection after pancreaticoduodenectomy. HBP (Oxford). 2016;18:893–899.

Read the full study here



Journal: Annals of Surgery

Title: Negative Pressure Wound Therapy for Surgical-site Infections: A Randomised Trial

Published: 1 June 2019

Javed AA, Teinor J, Wright M, Ding D, Burkhart RA, Hundt J, Cameron JL, Makary MA, He J, Eckhauser FE, Wolfgang CL, Weiss MJ. Negative Pressure Wound Therapy for Surgical-site Infections: A Randomised Trial. *Ann Surg.* 2019 Jun;269(6):1034–1040. doi: 10.1097/SLA.000000000000003056. **PMID: 31082899.**

Reduction of wound complication risk and length of stay with 3M™ Prevena™ Therapy.

Licari L, Campanella S, Carolla C, Viola S, Salamone G. Closed incision negative pressure therapy achieves better outcome than standard wound care: clinical outcome and cost-effectiveness analysis in open ventral hernia repair with synthetic mesh positioning. Cureus. 2020. 12(5):e8283.

Study design

Retrospective comparative cohort study (Italy).

Study purpose

The purpose of the study was to evaluate closed incision negative pressure therapy (ciNPT), Prevena Therapy, to standard of care (SOC) in regard to post-operative clinical outcomes and economical benefits for use in ventral hernia repair (VHR) with synthetic mesh positioning.

Methods

- ► Patients who underwent elective open VHR with synthetic mesh positioning from January 2015 to December 2017 at a single centre in Italy.
- ► Prevena Therapy (n=70) v. SOC (n=110).
- Patients followed for 90 days postoperatively.
- ► High Risk Inclusion Criteria: ≥ 1 risk factor.
 - ► Age > 65
 - ► Pre-existing wound infection
 - ► Pulmonary disease
 - ► BMI > 25
 - ► Malnutrition
 - ► Ascites
 - ► Hypertension

- ▶ Diabetes
- ► Active smoking
- ► Previous radiation therapy
- ► Steroid use
- ► Immunosuppression
- ► Chronic inflammatory disease

Summary of findings

The use of Prevena Therapy in high-risk populations following VHR with synthetic mesh significantly decreased the rate of complications and reduced the length of stay which resulted in a positive economic outcome.

Results Reduction in major Reduction in complications* superficial infections* 12 8^u (9/70) Prevena Therapy vs 4 3^u (3/70) Prevena Therapy vs. 43.6" (48/110) SOC (p<0.00001)* 22.7" (25/110) SOC (p=0.0006)* Reduction in mean Reduction in mean total cost per patient* in-hospital length of stay* 3 ± 1.37 Prevena Therapy vs. Prevena inpatient cost: € 4,230 6 ± 2.39 Control (p<0.00001)* SOC inpatient cost: 5,695 € (p=0.02)* Patient Cost Saving: € 1,465 Reduction in wound Reduction in deep infections* dehiscence¹ O^u (0/70) Prevena Therapy vs. 2.9^u (2/70) Prevena Therapy vs. 6.4" (7/110) SOC (p=0.04)*1 3.6^{u} (4/110) SOC (p=0.7)

Calculation(s) are derived based on relative patient group incidence rate reported in this study.

- *Statistically significant (p<0.05).
- †Note: The use of Prevena Therapy for the reduction in the incidence of deep infections and wound dehiscence has not been reviewed by the U.S. FDA.

The effectiveness of Prevena Therapy in reducing the incidence of SSIs and seroma in all surgical procedures and populations has not been demonstrated. See full indications for use and limitations at hcbgregulatory.3M.com.

Read the full study here



Journal: Cureus

Title: Closed Incision Negative Pressure Therapy Achieves Better Outcome Than Standard Wound Care: Clinical Outcome and Cost-Effectiveness Analysis in Open Ventral Hernia Repair With Synthetic Mesh Positioning

Published: 26 May 2020

Licari L, Campanella S, Carolla C, Viola S, Salamone G. Closed incision negative pressure therapy achieves better outcome than standard wound care: clinical outcome and cost-effectiveness analysis in open ventral hernia repair with synthetic mesh positioning. Cureus. 2020. 12(5):e8283.

OPEN ACCESS

Reduction of the incidence of surgical site infection with 3M™ Prevena™ Therapy in emergency laparotomy patients.

Cheong Chung JN, Ali O, Hawthornthwaite E, Watkinson T, Blyth U, McKigney N, Harji DP, Griffiths B. Closed incision negative pressure wound therapy is associated with reduced surgical site infection after emergency laparotomy: A propensity matched-cohort analysis. *Surgery*. 2021 Nov;170(5):1568–1573.

Study design

Retrospective comparative cohort study (United Kingdom).

Study purpose

To evaluate with a propensity matched analysis whether the use of closed incision negative pressure therapy (ciNPT), Prevena Therapy, decreases surgical site infections (SSI) compared to standard surgical dressings after emergency laparotomy.

Methods

- A registry-based, cohort study was undertaken using data the NELA registry.
- ► The National Emergency Laparotomy Audit (NELA) is part of the National Clinical Audit and Patient Outcomes Programme (NCAPOP), overseen by the Healthcare Quality Improvement Partnership (HQIP) in the UK.
- ▶ 1484 patients identified from the NELA dataset.
- Propensity score matching resulted in two equally matched cohorts with 237 patients in each arm.
- ► Prevena Therapy applied of midline incision and left in situ for 7 days or until discharge if before.
- ► Standard surgical dressing (Opsite dressing).
- Primary outcome was SSI per Centers for Disease Control criteria.

Summary of findings

- ► This registry-based cohort study using the NELA registry uses real world data to shows the use of Prevena Therapy in emergency laparotomy patients is associated with a significant reduction of surgical site infections (33.8% vs 16.9%; p<0.001*).
- ► The study also demonstrated a reduction in both superficial and deep SSI.

Calculation(s) are derived based on relative patient group incidence rate reported in this study. *Statistically significant (p<0.05).

†Note: The use of Prevena Therapy for the reduction in the incidence of deep infections and wound dehiscence has not been reviewed by the U.S. FDA.



(continued)

(continued)

Read the full study here



Journal: Surgery

Title: Closed incision negative pressure wound therapy is associated with reduced surgical site infection after emergency laparotomy: a propensity matched-cohort analysis.

Published: 26 May 2021

Cheong Chung JN, Ali O, Hawthornthwaite E, Watkinson T, Blyth U, McKigney N, Harji DP, Griffiths B. Closed incision negative pressure wound therapy is associated with reduced surgical site infection after emergency laparotomy: A propensity matched-cohort analysis. Surgery. 2021 Nov;170(5):1568-


1573. OPEN ACCESS

PRM in general surgery Decision guide

Decision guide

Patient and procedure risk stratification in gneral surgery backed by clinical evidence.


While most surgical patients may benefit from 3M™ Prevena™ Therapy, patients at high risk for complications such as surgical site infection may see added benefit. The following uses select study data¹-8 to provide an illustrative guide to aid in risk stratification. This is not an all-inclusive list of risk factors. Clinicians are advised to use their clinical judgment to identify high-risk patients or high-risk procedures.



References

1. Willy C, Agarwal A, Andersen CA, De Santis G, Gabriel A, Grauhan O, Guerra OM, Lipsky BA, Malas MB, Mathiesen LL, Singh DP, Reddy VS. Closed incision negative pressure therapy: international multidisciplinary consensus recommendations. *Int Wound J.* 2017 Apr;14(2):385-398. **OPEN ACCESS** 2. Curran T, Alvarez D, Pastrana Del Valle J, Cataldo TE, Poylin V, Nagle D. Prophylactic closed incision negative pressure wound therapy is associated with decreased surgical site infection in high-risk colorectal surgery laparotomy wounds. *Colorectal Dis.* 2019 Jan. 21(1):110-118. **OPEN ACCESS** 3. Javed AA, Teinor J, Wright M, Ding D, Burkhart RA, Hundt J, Cameron JL, Makary MA, He J, Eckhauser FE, Wolfgang CL, Weiss MJ. Negative Pressure Wound Therapy for Surgical-site Infections: A Randomised Trial. *Annals of Surgery.* 2019 Jun;269(6):1034-1040. **PMID 31082899** 4. Zaidi A, El-Masry S. Closed incision negative pressure therapy in high-risk general surgery patients following laparotomy: a retrospective study. *Colorectal Disease.* 2017 Mar;19(3):283-287. **OPEN ACCESS** 5. Licari L, Campanella S,

Carolla C, Viola S, Salamone G. Closed Incision Negative Pressure Therapy Achieves Better Outcome Than Standard Wound Care: Clinical Outcome and Cost-Effectiveness Analysis in Open Ventral Hernia Repair With Synthetic Mesh Positioning. Cureus. 2020;12(5):e8283. OPEN ACCESS 6. Ayuso SA, Elhage SA, Okorji LM, Kercher KW, Colavita PD, Heniford BT, Augenstein VA. Closed-Incision Negative Pressure Therapy Decreases Wound Morbidity in Open Abdominal Wall Reconstruction With Concomitant Panniculectomy. Ann Plast Surg. 2022 Apr 1;88(4):429-433. PMID 34670966 7. Cheong Chung JN, Ali O, Hawthornthwaite E, Watkinson T, Blyth U, McKigney N, Harji DP, Griffiths B. Closed incision negative pressure wound therapy is associated with reduced surgical site infection after emergency laparotomy: A propensity matched-cohort analysis. Surgery. 2021 May 26:S0039-6060(21)00334-2. PMID 34052025 8. Lakhani A, Jamel W, Riddiough GE, Cabalag CS, Stevens S, Liu DS. Prophylactic negative pressure wound dressings reduces wound complications following emergency laparotomies: A systematic review and meta-analysis. Surgery. 2022 Sep;172(3):949-954. PMID 35779950



3M[™] Prevena[™] Therapy Units



3M™ Prevena™ 125 Therapy Unit



3M™ Prevena™ Plus 125 Therapy Unit (7 or 14 day)

3M™ Prevena™ Dressings

3M™ Prevena™ Peel and Place Dressings



3M™ Prevena™ Peel and Place Dressing -13 cm



3M™ Prevena™ Peel and Place Dressing - 20 cm



3M™ Prevena™ Peel and Place Dressing – 35 cm

3M™ Prevena™ Plus **Customizable Dressing**



3M™ Prevena™ Plus Customizable Dressing

3M™ Prevena Restor™ Dressings



3M™ Prevena Restor™ Arthro·Form™ Dressing



3M™ Prevena Restor™ Bella·Form™ Dressing

Note: The availability of Prevena Therapy Units and Dressings may vary by country. To get precise information about availability, we recommend reaching out to your local sales representative.



Healthcare professionals:

Besuchen Sie <u>3M.de/Prevena</u> um mehr zu erfahren oder eine Produktberatung anzufordern.

Hinweis: Für das Prevena Incision Management System gelten spezifische Indikationen, Kontraindikationen, Warnungen, Vorsichtsmaßnahmen und Sicherheitsinformationen. Bitte machen Sie sich vor der Applikation mit der Bedienungsanleitung für Anwender des Prevena Systems vertraut. Dieses Material ist für medizinisches Fachpersonal bestimmt.

3M[™] Prevena Therapie und 3M[™] Prevena Restor Therapie Gebrauchsindikation:

Das Prevena Incision Management System und das Prevena Restor Incision Management System sind für das Management der Umgebung geschlossener chirurgischer Inzisionen und der umgebenden intakten Haut bei Patienten bestimmt, bei denen ein Risiko für postoperative Komplikationen, wie z. B. einer Infektion, besteht. Dabei wird eine geschlossene Umgebung durch die Anwendung eines Unterdruck-Wundtherapiesystems auf der Inzision aufrechterhalten. Die Prevena Dressings mit Silber reduzieren die mikrobielle Besiedlung des Gewebes.

* Die Wirksamkeit der Prevena Therapie bei der Reduzierung der Inzidenz von postoperativen Wundinfektionen und Seromen wurde nicht bei allen chirurgischen Eingriffen und Patientengruppen nachgewiesen. Die vollständigen Indikationen zur Anwendung und Einschränkungen finden Sie auf <u>hcbgregulatory.3M.com</u>

Dieses Dokument wurde von 3M erstellt. Es umfasst nicht alle ciNPT-Veröffentlichungen.





3M Deutschland GmbH Health Care Business Carl-Schurz-Str. 1 41453 Neuss

Tel +49 (0)2131 881 9260 3Mmedica.de@mmm.com www.3M.de/medicalsolutions **3M (Schweiz) GmbH Health Care** Eggstrasse91 8803 Rüschlikon

Tel +41 (0)44 724 90 90 innovation.ch@mmm.com www.3Mschweiz.ch/medicalsolutions

3M Österreich GmbH Health Care Kranichberggasse 4 1120 Wien

Tel +43 (0)186 68 60 medizin-at@mmm.com www.3Maustria.at/medicalsolutions

© 3M 2023. 3M, Prevena und Prevena Restor sind eingetragene Marken oder Warenzeichen der 3M Company. Alle Rechte vorbehalten. d6d0312760af