solventum
V.A.C.® Therapy

Reminders for V.A.C.® Therapy with V.A.C.® Peel and Place Dressing

Indications

V.A.C.® Peel and Place Dressing can be used on:

- Chronic wounds
- Acute wounds
- Traumatic wounds
- Subacute wounds
- Dehisced wounds
- Partial-thickness burns
- Ulcers (such as diabetic, pressure or venous insufficiency)
- Flaps
- Grafts

Warnings

Do not use on wounds with:



- Undermining greater than 2 cm
- Tunneling
- Depth greater than 6 cm

Dressing should cover the entire wound bed and extend onto intact skin.

Find the right size

Small - EZ5SML	Medium - EZ5MED	Large - EZ5LRG
Max wound size: 6.1 cm x 8.6 cm **Max depth:** 2 cm	**Max wound size:** 11.1 cm x 16.6 cm **Max depth:** 4 cm	**Max wound size:** 13.6 cm x 24.2 cm **Max depth:** 6 cm

Always read and follow the detailed Instructions for Use along with important safety information provided in the packaging.

solventum
V.A.C.® Therapy

Application tips for V.A.C.® Therapy with V.A.C.® Peel and Place Dressing

Wear time

Dressing can remain in place for **up to 7 days**.[1]

Device settings

Set the device to **-75 to -150 mmHg** on **continuous mode**.

Do Not	Drape Tips
• Cut foam. • If skin barrier/protectant is used under the adhesive portion of the dressing, do not reposition the dressing as adhesive properties may be reduced. • Use with other V.A.C.® Therapy Dressing Kits or foams, as they have a shorter wear time. • Push foam portion of the dressing down. Allow the dressing to draw down and conform the wound bed as negative pressure is applied. • Stretch dressing.	• Apply loosely over the wound area. • Leave at least a 5 cm border of drape wherever possible. • Can be lifted and repositioned once during initial application. • For curved anatomy, slits may be cut in the drape portion of the dressing to help reduce overlap, drape folds or remove wrinkles.

Solventum Malaysia
Level 8A, Block F, Oasis Square,
No. 2, Jalan PJU 1A/7A, Ara Damansara,
47310 Petaling Jaya, Selangor, Malaysia.

Watch how-to-apply video →

Note: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product Instructions for Use prior to application. This material is intended for healthcare professionals.

[1] SAT-BSER-05-869347 VAC Peel and Place BSER. 510(k) K222859.

This material is only meant for Healthcare Professionals. Registered under Act 737. MDA Registration number: GC7625226-235928