



The Bowie-Dick revolution begins

3M™ Attest™ eBowie-Dick Test System

When the pressure is on, stop guessing and start knowing

A lot is dependent on your steam sterilizer. You need pass/fail results quickly and definitively, without guesswork, while improving workflow. The 3M™ Attest™ eBowie-Dick Test System revolutionizes sterilization testing and delivers automated results. The same testing method you've always trusted, now digital.



No more guesswork

Provides quick, definitive pass/fail results of your steam sterilizer, removes subjective interpretations and provides more assurance and confidence.



Enhances workflow

Removes the burden of paperwork and eliminates manual recordkeeping — save time and operate more efficiently.



Innovative technology

First-of-its-kind system backed by Solventum science — so innovative, it's in a category of its own created by the FDA.

Simplify. Standardize. Streamline.



Designed for sustainability

The 3M™ Attest™ eBowie-Dick Test System is part of our Strategic Sustainability Framework, which can help drive your own sustainability efforts.

- ✓ Recyclable card holder is made of Polyetherimide plastic.*
- ✓ Compliant with global electronic waste directives RoHS and WEEE.
- ✓ eBowie-Dick test cards can reduce waste entering a hospital's waste stream by 94%.†
- ✓ No intentionally added natural rubber latex in the product or primary packaging.
- ✓ Can eliminate waste removal need for approximately 88 lbs. of paper waste per sterilizer per year.‡

*Recycling programs for this product may not exist in your area.

†Compared to waste generated by 3M Bowie-Dick paper test packs. Waste mass based on eBowie-Dick 10135B product and packaging case weight, with percent calculated on a per test basis.

‡Per AAMI ST79 usage recommendations and 3M estimate of reduction of waste when adopting eBowie-Dick.

Component	Catalog no.	Description	Packaging	Shipping
3M™ Attest™ eBowie-Dick Test Card	10135	132°C (270°F), 134°C (273°F), and 135°C (275°F)	30 cards/bag	1 bag/case
3M™ Attest™ eBowie-Dick Test Card — Bulk	10135B	132°C (270°F), 134°C (273°F), and 135°C (275°F)	30 cards/bag	4 bags/case
3M™ Attest™ eBowie-Dick Card Holder	10135CH		2 card holders/box	1 box/case
3M™ Attest™ eBowie-Dick Card Holder — Bulk	10135CHB		6 card holders/box	1 box/case
3M™ Attest™ eBowie-Dick Auto-reader	1190	12,9cm (W) x 6,3cm (H) x 4,7cm (D) [5.8" (W) x 2.5" (H) x 1.9" (D)]	1 unit/box 2 card holders/box	1 box/case

Available in Canada from:



3M Healthcare Canada Company Limited
75 Courtneypark Dr W, Unit #4
Mississauga
ON L5W 0E3
Canada
Phone: 1-833-954-0525
Web: 3M.ca/Medical

Contact your Solventum representative

Everyday support: Solventum Medical Solutions Service Line
1-833-954-0525

The indications are provided in the IFU for the test card.

The indications for use should align with the sterilization cycle to be monitored.

© Solventum 2024. Solventum, the S logo and Attest are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. Other trademarks are the property of their respective owners. CA-EN-70-2013-1655-4

Learn more

