




No assembly required

3M™ Attest™ Super Rapid Vaporized
Hydrogen Peroxide Clear Challenge Pack



Elevate your VH2O2 sterilization monitoring practice

Save time and reduce potential error with the 3M™ Attest™ Super Rapid Vaporized Hydrogen Peroxide Clear Challenge Pack including both a BI and a CI. It's the first and only preassembled VH2O2 test pack that can be used for routine monitoring and qualification across multiple sterilizer brands, models and cycle types.*

3M Health Care is now Solventum

Streamline sterilization workflow and inventory



The first ever preassembled VH2O2 test pack contains both a CI and BI

- ✓ Saves time and reduces chance for error—no need to assemble peel pouches
- ✓ Elevates your VH2O2 sterilization monitoring workflow to a similar level as steam sterilization best practices
- ✓ Preassembled design more accurately simulates the environment inside a surgical pack than using a stand-alone BI for load monitoring
- ✓ Validated across multiple sterilizer brands, models and cycle types, helping you save on inventory

Ordering information

Product	Catalog no.	Packaging	Shipping
3M™ Attest™ Super Rapid Vaporized Hydrogen Peroxide Clear Challenge Pack	1295PCD-24	6 PCDs/bag	4 bags and 24 controls/case
3M™ Attest™ Super Rapid Vaporized Hydrogen Peroxide Clear Challenge Pack	1295PCD-05	6 PCDs/bag	4 bags and 5 controls/case



Talk to your local Solventum Representative to find out more.



Solventum Medical Surgical
10 Ang Mo Kio Street 65
Techpoint #01-01
Singapore 569059

Phone +65 6577 1266
Email: sales.sg@solventum.com

*as of November 11, 2024
Solventum, formerly 3M Health Care
©Solventum 2025. Solventum, the S logo and Attest are trademarks of Solventum or its affiliates. 3M, and the 3M logo are trademarks of 3M. Other trademarks are the property of their respective owners. US-70-2013-1753-7