

Decision guide

3M™ Ioban™ 2 Antimicrobial Incise Drape and
3M™ Ioban™ CHG Chlorhexidine Gluconate
Incise Drape (contains 2% w/w CHG)



Use these charts to see which
antimicrobial incise drape may
be appropriate for your surgical
application based on product
features and dimensions



Orthopedic surgical procedures

Use this chart to see which antimicrobial incise drape may be appropriate for your surgical application based on product features and dimensions.

3M Cat. No.	Product	Adhesive area	Items/box	Boxes/case	ACL	Arthroscopy	ORIF small extremity	Hip pinning (prox. or distal)	Hip pinning (interlocking nail)	Pediatrics	Podiatry (including ankle)	Spine	Total hip	Total hip (anterior)	Total knee	Total shoulder
	6635	3M™ Ioban™ 2 Antimicrobial Incise Drape	3.87" x 7.875" 10cm x 20cm	10	4											
	6640	3M™ Ioban™ 2 Antimicrobial Incise Drape	13" x 13" 34cm x 35cm	10	4											
	6640EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	13" x 13" 35cm x 35cm	10	4											
	8840EZ	3M™ Ioban™ CHG Chlorohexidine Gluconate Incise Drape	13" x 13" 34cm x 35cm	10	4											
	6650	3M™ Ioban™ 2 Antimicrobial Incise Drape	22" x 17" 56cm x 45cm	10	4											
	6650EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	23" x 17" 60cm x 45cm	10	4											
	8850EZ	3M™ Ioban™ CHG Chlorohexidine Gluconate Incise Drape	22" x 17" 56cm x 45cm	10	4											
	6648	3M™ Ioban™ 2 Antimicrobial Incise Drape	22" x 23" 56cm x 60cm	10	4											
	6648EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	23" x 23" 60cm x 60cm	10	4											
	6651	3M™ Ioban™ 2 Antimicrobial Incise Drape	22" x 33" 56cm x 85cm	10	4											
	6651EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	23" x 33" 60cm x 85cm	10	4											
	8851EZ	3M™ Ioban™ CHG Chlorohexidine Gluconate Incise Drape	22" x 33" 56cm x 85cm	10	4											
	6617	Isolation Drape with 3M™ Ioban™ 2 Antimicrobial Incise Film and Pouch	Overall size: 125" x 83" / 320cm x 213cm Adhesive size: 19" x 9.37" / 50cm x 24cm	5	4											
	6619	Large Isolation Drape with 3M™ Ioban™ 2 Antimicrobial Incise Film and Pouch	Overall size: 129" x 100" / 328cm x 254cm Adhesive size: 27" x 12" / 70cm x 32cm	5	1											

 Incise drape for your application  Secondary option  Specialty incise drape

OB/GYN surgical procedures

Use this chart to see which antimicrobial incise drape may be appropriate for your surgical application based on product features and dimensions.

	3M Cat. No.	Product	Adhesive area	Items/box	Boxes/case	Open abdominal/ pelvic surgery	C-section
	6650	3M™ Ioban™ 2 Antimicrobial Incise Drape	22" x 17" 56cm x 45cm	10	4		
	6650EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	23" x 17" 60cm x 45cm	10	4		
	8850EZ	3M™ Ioban™ CHG Chlorohexidine Gluconate Incise Drape	22" x 17" 56cm x 45cm	10	4		
	6648	3M™ Ioban™ 2 Antimicrobial Incise Drape	22" x 23" 56cm x 60cm	10	4		
	6648EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	23" x 23" 60cm x 60cm	10	4		
	8851EZ	3M™ Ioban™ CHG Chlorohexidine Gluconate Incise Drape	22" x 33" 56cm x 85cm	10	4		
	6657	3M™ Steri-Drape™ Pouch with Ioban™ 2 Incise Film	Overall size: 34" x 29" / 89cm x 76cm Incise area: 11" x 11" / 30cm x 30cm	10	4		
	6658	3M™ Steri-Drape™ Pouch with Ioban™ 2 Incise Film	Overall size: 30" x 30" / 76cm x 76cm Incise area: 13" x 17" / 33cm x 43cm	5	4		
	6659	3M™ Steri-Drape™ Pouch with Ioban™ 2 Incise Film	Overall size: 29" x 34" / 74cm x 87cm Incise area: 16" x 20" / 43cm x 52cm	5	4		

 Incise drape for your application  Secondary option

General surgical procedures

Use this chart to see which antimicrobial incise drape may be appropriate for your surgical application based on product features and dimensions.

3M Cat. No.	Product	Adhesive area	Items/box	Boxes/case	Abdominal-perineal resection	Open appendectomy	Open colon resection	Open hernia repair	Laparotomy	Liver transplant	Kidney transplant
	6661EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	10.5" x 8" 26cm x 20cm	50	2						
	6640	3M™ Ioban™ 2 Antimicrobial Incise Drape	13" x 13" 34cm x 35cm	10	4						
	6640EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	13" x 13" 35cm x 35cm	10	4						
	8840EZ	3M™ Ioban™ CHG Chlorohexidine Gluconate Incise Drape	13" x 13" 34cm x 35cm	10	4						
	6650	3M™ Ioban™ 2 Antimicrobial Incise Drape	22" x 17" 56cm x 45cm	10	4						
	6650EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	23" x 17" 60cm x 45cm	10	4						
	8850EZ	3M™ Ioban™ CHG Chlorohexidine Gluconate Incise Drape	22" x 17" 56cm x 45cm	10	4						
	6648	3M™ Ioban™ 2 Antimicrobial Incise Drape	22" x 23" 56cm x 60cm	10	4						
	6648EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	23" x 33" 60cm x 60cm	10	4						

 Incise drape for your application  Secondary option

Vascular/cardiovascular/thoracic surgical procedures

Use this chart to see which antimicrobial incise drape may be appropriate for your surgical application based on product features and dimensions.

3M Cat. No.	Product	Adhesive area	Items/box	Boxes/case	Open abdominal aortic aneurysm (AAA)	AV fistula	Coronary artery bypass (CAB)	Coronary artery bypass with graft (CABG)	Carotid endarterectomy	Embolectomy	Femoropopliteal bypass	Lobectomy	Valve replacement	Pacemaker
	6635	3M™ Ioban™ 2 Antimicrobial Incise Drape	3.875" x 7.875" 10cm x 20cm	10	4									
	6640	3M™ Ioban™ 2 Antimicrobial Incise Drape	13" x 13" 34cm x 35cm	10	4									
	6640EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	13" x 13" 35cm x 35cm	10	4									
	8840EZ	3M™ Ioban™ CHG Chlorohexidine Gluconate Incise Drape	13" x 13" 34cm x 35cm	10	4									
	6661EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	10.5" x 8" 26cm x 20cm	50	4									
	6650	3M™ Ioban™ 2 Antimicrobial Incise Drape	22" x 17" 56cm x 45cm	10	4									
	6650EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	23" x 17" 60cm x 45cm	10	4									
	8850EZ	3M™ Ioban™ CHG Chlorohexidine Gluconate Incise Drape	22" x 17" 56cm x 45cm	10	4									
	6648	3M™ Ioban™ 2 Antimicrobial Incise Drape	22" x 23" 56cm x 60cm	10	4									
	6648EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	23" x 23" 60cm x 60cm	10	4									
	6651	3M™ Ioban™ 2 Antimicrobial Incise Drape	22" x 33" 56cm x 85cm	10	4									
	6651EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	23" x 33" 60cm x 85cm	10	4									
	8851EZ	3M™ Ioban™ CHG Chlorohexidine Gluconate Incise Drape	22" x 33" 56cm x 85cm	10	4									



Incise drape for your application



Secondary option

Neurologic/spinal surgical procedures

Use this chart to see which antimicrobial incise drape may be appropriate for your surgical application based on product features and dimensions.

3M Cat. No.	Product	Adhesive area	Items/box	Boxes/case	Craniotomy	VP shunt	Spine
 6640	3M™ Ioban™ 2 Antimicrobial Incise Drape	13" x 13" 34cm x 35cm	10	4			
 6640EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	13" x 13" 35cm x 35cm	10	4			
 8840EZ	3M™ Ioban™ CHG Chlorohexidine Gluconate Incise Drape	13" x 13" 34cm x 35cm	10	4			
 6650	3M™ Ioban™ 2 Antimicrobial Incise Drape	22" x 17" 56cm x 45cm	10	4			
 6650EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	23" x 17" 60cm x 45cm	10	4			
 8850EZ	3M™ Ioban™ CHG Chlorohexidine Gluconate Incise Drape	22" x 17" 56cm x 45cm	10	4			
 6648	3M™ Ioban™ 2 Antimicrobial Incise Drape	22" x 23" 56cm x 60cm	10	4			
 6648EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	23" x 23" 60cm x 60cm	10	4			
 6651	3M™ Ioban™ 2 Antimicrobial Incise Drape	22" x 33" 56cm x 85cm	10	4			
 6651EZ	3M™ Ioban™ 2 Antimicrobial Incise Drape EZ	23" x 33" 60cm x 85cm	10	4			
 8851EZ	3M™ Ioban™ CHG Chlorohexidine Gluconate Incise Drape	22" x 33" 56cm x 85cm	10	4			

 Incise drape for your application  Secondary option

For more information about how our antimicrobial incise drapes can help reduce the risk of surgical site infections, contact your Solventum account representative or visit go.solventum.com/lobanPDF.



2510 Conway Ave.
St. Paul, MN 55144 U.S.A

Phone 800 228-3957
Web Solventum.com

© Solventum 2025. Solventum and the S logo are trademarks of Solventum or its affiliates. 3M is a trademark of 3M company. Other trademarks are the property of their respective owners.