




Remove uncertainty with Every Load Monitoring (ELM)


Trusted partner in patient safety

ANSI/AAMI ST79 Standards set the baseline of the use of a biological indicator (BI) for routine monitoring at least weekly, preferably daily, and every implant load. For Vaporized Hydrogen Peroxide, ANSI/AAMI ST58 standards state a BI should be used daily, preferably every load, and every implant load.

We'd love to be your ELM partner so **every** load is monitored with a BI and **every** patient gets the highest level of sterility assurance each day.



Monitoring all steam **or** vaporized hydrogen peroxide (VH2O2) loads



Monitoring all steam **and** VH2O2 loads



How does ELM work?

Commit to two years of Solventum sterilization assurance product purchases* and we will send you a custom engraved ELM plaque to display your patient safety dedication. For each year of participation, we'll send annual magnets to add to your plaque.

*Requires you to agree to a two-year commitment to purchase any combination of biological indicators and test packs for steam and/or vaporized hydrogen peroxide (VH2O2) to support your ELM practice. The commitment will end two years thereafter unless otherwise terminated by your facility upon ninety (90) days written notice to Solventum. Early termination will require your facility to return the customized engraved plaque to Solventum. This offer is available to health care provider facilities with sterilization departments within the 50 states that comprise the United States of America and Washington D.C.

What products are eligible?

Any combination of biological indicators and test packs for steam and/or vaporized hydrogen peroxide (VH2O2) products, to support your ELM practice, during the two year period.



Method	Catalog no.	Product	Description / Indications for use	Packaging	Shipping
 Steam	1492PCD-24	3M™ Attest™ Super Rapid Steam Clear Challenge Pack	Dynamic-air-removal (pre-vacuum and SFPP): 132°C (270°F) = 4 minutes 134°C (273°F) = 4 minutes 135°C (275°F) = 3 minutes	24 packs	24 controls/case
	1492PCD-05			24 packs	5 controls/case
	1493PCDG-24	3M™ Attest™ Super Rapid Steam Gravity Clear Challenge Pack	Gravity displacement: 121°C (250°F) = 30 minutes	24 packs	24 controls/case
	1492PCDE-24	3M™ Attest™ Super Rapid Steam Extended Cycle Clear Challenge Pack	Dynamic-air-removal: 132°C (270°F) = 10 minutes	24 packs	24 controls/case
	1493	3M™ Attest™ Super Rapid Readout Biological Indicator	Gravity: 121°C (250°F) = 30 minutes Dynamic-air-removal (pre-vacuum and SFPP): 121°C (250°F) = 15, 20, 30 or 35 minutes	50 each/bag	4 bags/case
	1493P			20 each/bag	5 bags/case
	1491	3M™ Attest™ Super Rapid Readout Biological Indicator	Gravity: 132°C (270°F) = 3 or 10 minutes 135°C (275°F) = 3 or 10 minutes	50 each/bag	4 bags/case
	1492V	3M™ Attest™ Super Rapid Readout Biological Indicator	Dynamic-air-removal (pre-vacuum and SFPP): 132°C (270°F) = 3, 3.5, 4, 5.5 or 6 minutes 134°C (273°F) = 3 or 4 minutes 135°C (275°F) = 3, 3.5 or 10 minutes	50 each/bag	4 bags/case
 VH2O2	1295PCD-24	3M™ Attest™ Super Rapid Vaporized Hydrogen Peroxide Clear Challenge Pack	Vaporized Hydrogen Peroxide (VH2O2)*	6 packs/bag	4 bags and 24 controls/case
	1295PCD-05			6 packs /bag	4 bags and 5 controls/case
	1295	3M™ Attest™ Rapid Readout Biological Indicator	Vaporized Hydrogen Peroxide (VH2O2)*	30 each/bag	4 bags/case

*Refer to the indications for use provided in the biological indicator or PCD Instructions for Use. The indications for use should align with the sterilization cycle to be monitored.

The value of the customized engraved plaque that you are receiving from Solventum is \$109.84 under federal law at 42 U.S.C. §1320a-7b(b)(3) (A) and regulations issued there under. Accordingly, you are responsible for reporting/or providing information on all discounts, rebates, or reductions in price to reimbursement agencies (including Medicare and Medicaid) and other entities in accordance with all applicable laws and regulations, including the provisions of the discount exception at 42 U.S.C. §1320a-7b(b)(3)(A) and the discount safe harbor regulation at 42 C.F.R. §1001.952(h). Further, you shall promptly disclose and appropriately reflect any discount, rebate or other reduction in price as may be required under Medicare or any other federal or state health program, consistent with the requirements of the Medicare/Medicaid Antikickback Statute and regulations.



Solventum Corporation
1750 Yankee Doodle Road
Eagan, MN 55121

Phone 800-228-3975
Web Solventum.com

Contact your Solventum representative

Everyday support: Solventum Medical Service Line | 1-800-228-3957/Attest

Solventum, formerly 3M Health Care

© Solventum 2026. Solventum, the S logo and Attest are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. 70-2013-1994-7