



ECG Electrodes product comparison chart

Red Dot™ Electrodes

The Right Electrode Matters

Choosing the right electrode can be difficult. That's why we've optimised our 3M™ Red Dot™ ECG Electrode portfolio to feature the options that meet your most common patient and performance needs. With 3M, optimal adhesion for a variety of patients, consistent and reliable trace quality, and ease of use come standard—helping you make a positive impact on patient care. **Reach for the Red.●**



	Adult				Adult/Pediatric		Neonatal	
	Universal				Diagnostic			
	Diaphoretic	Diaphoretic	Repositionable		Gentle			
	2228 2244	2560 2570	2660 2670	2239	2268	2330 2360*	2269T	2282
Adhesion			2660: 2670:			2330: 2360:		
Trace Quality								
MR Conditional	2228: No 2244: Yes	2560: No 2570: Yes	Yes	No	Yes	Yes	Yes	No
Radiolucent	2228: No 2244: Yes	2560: No 2570: Yes	Yes	No	Yes	Yes	Yes	Yes
Backing	Foam	Foam	Soft Cloth	Micropore Tape	Soft Cloth	Polyester Film	Clear Tape	Soft Cloth
Pack Size	50	50	3, 5	50	3	100	3	3

*Adult only

3M United Kingdom PLC
Charnwood Campus
10 Bakewell Road
Loughborough
LE11 5RB
01509 611611
3M.co.uk/RedDot

3M Ireland Limited
The Iveagh Building
The Park, Carrickmines
D18 X015
Ireland
+353 (0)1 280 3555
3M.ie/RedDot

3M and Red Dot are trademarks of 3M. © 3M 2018. All Rights Reserved.