



A simplified solution worth looking into

3M™ Attest™ Steam Clear Challenge Packs



Meet the most comprehensive solution for steam sterilization monitoring*

Gain confidence in your steam sterilization cycles with our new, easy-to-use clear challenge packs. They're preassembled and allow you to inspect the BI and the CI before and after processing without opening the pack.

*3M™ Attest™ Steam Clear Challenge Packs (1492PCD, 1492PCDE, 1493PCDG) offer customers 8 cycle claims cleared by US FDA across the 3 challenge packs combined. Based on U.S. FDA cycle clearances as of January 17, 2025.

Clearly a better way to monitor sterilization cycles

3M™ Attest™ Steam Clear Challenge Packs



Preassembled and easy to use





Transparent design allows you to inspect the BI and CI before and after processing without opening the pack



Represents the challenge posed by instruments sterilized every day

May reduce Process Challenge Device (PCD) waste up to 82% by weight compared to the 3M™ Attest™ Super Rapid 5 Steam-Plus Challenge Pack*

Steam ordering information

STEAM	Product	Catalogue no.	Packaging	Shipping
	3M™ Attest™ Super Rapid Steam Clear Challenge Pack	1492PCD-24	24 packs	24 controls/case
	3M™ Attest™ Super Rapid Steam Clear Challenge Pack	1492PCD-05	24 packs	5 controls/case
	3M™ Attest™ Super Rapid Steam Gravity Clear Challenge Pack	1493PCDG-24	24 packs	24 controls/case
	3M™ Attest™ Super Rapid Steam Extended Cycle Clear Challenge Pack	1492PCDE-24	24 packs	24 controls/case

	Biological Indicators	Catalogue no.	Packaging	Shipping
NEW	3M™ Attest™ Super Rapid Readout Biological Indicator	1493	50 each/bag	4 bags/case
	3M™ Attest™ Super Rapid Readout Biological Indicator	1493P	20 each/bag	5 bags/case
	3M™ Attest™ Super Rapid Readout Biological Indicator	1491	50 each/box	4 boxes/case
	3M™ Attest™ Super Rapid Readout Biological Indicator	1492V	50 each/box	4 boxes/case

*"Waste" refers only to the weight of discarded product and packaging. This is not an environmental impact claim. Comparison assumes equivalent use per sterilization cycle. Other factors like recyclability, emissions or lifecycle performance are not included. Supporting data available upon request.



Solventum Medical Surgical
75 Courtneypark Drive West
Unit #4, Mississauga, ON
L5W 0E3

Contact your Solventum representative

Everyday support: Solventum Medical Service Line

1-800-222-3957

Solventum, formerly 3M Health Care

© Solventum, 2025. Solventum, the S logo and Attest are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. Other trademarks are the property of their respective owners.

70-2013-1770-1-EN-CA